

ZEBRA INVENTO™

Positioning chair

ZEBRA Invento™ is a positioning chair designed for children and adolescents who need support in the process of sitting independently or in adopting and maintaining a sitting position. The device combines the function of a seat and a lying position, which combined with the smooth adjustment of tilt-in-space function, supports the rehabilitation process. Adjustable headrest, footrests, chest and pelvic supports help maintain a correct lying, semi-lying or sitting position. The functionality of ZEBRA Invento™ makes the positioning chair meet all the expectations of even the most demanding users, significantly increasing everyday comfort, convenience and safety.



Tilt-in-space
function

24⁺

24-month
warranty




Stable
construction



Easy
operation

Colors


ATMO™ upholstery
Breathable and comfortable



Grass



Red



Sea



Grey



The tilt-in-space function can be easily achieved by the adjustment of the frame.



Lying position – a unique function of the seat ensuring full comfort. With an extended backrest and/or deeper seat, additional filling pillows guarantee comfortable back support.




70 - 150 cm

Indicative user
height



25 - 60 kg

Max. user
weight



For **more details** about
the product, please visit
our website.

This is a medical device.
Use it in accordance with the user manual or the label.

INDICATIONS

Cerebral palsy (CP) – various types and forms | Brain and spinal cord injuries of various etiologies with limb paresis, severe balance disturbances, and coordination impairments | Demyelinating diseases (multiple sclerosis – MS, amyotrophic lateral sclerosis – ALS, and others) | Meningocele | Muscle diseases with paresis and paralysis (including dystrophy, myasthenia gravis, polymyositis) | Genetic syndromes with limb paresis | Degenerative and metabolic diseases of the nervous system | Degenerative and rheumatic joint diseases | Other diseases with limb paresis, paralysis, and musculoskeletal impairments | After injuries, following traumatic brain injuries or spinal cord injuries with limb paresis and paralysis, particularly with quadriplegia (as part of the rehabilitation process or as permanent orthopedic support)

CONTRAINDICATIONS

Observation: Monitor user during standing. Stop therapy if paleness, sweating, fainting, joint overload, or body deformations occur; consult with a doctor or physiotherapist | **Absolute Contraindication:** Critical condition of the user. | **Relative Contraindications:** Resting hypertension > 200/110 mmHg, pulmonary hypertension, heart rhythm disorders, unstable asthma, unstable metabolic issues, moderate heart valve disease. | Conditions like healing fractures or severe osteoporosis that prevent weight-bearing | Incomplete fracture healing (unless approved by a doctor) | Postoperative joint load restrictions | Pain from dislocation, excessive tissue stretch, or pressure; dislocation without pain is not contraindicated | Significant hip/knee flexion or ankle contractures; ensure no excessive pressure or stretch on affected areas.