

# OMBRELO™

Special needs stroller

OMBRELO™ is a lightweight, easy-to-fold special needs stroller for children and adults. It features the innovative BodyPrint™ seat with a unique belt system that ensures optimal positioning and comfort, even for extended use. The adjustable seat depth and backrest height provide a tailored fit, while the 27° fixed tilt aids upper body positioning for those with limited trunk stability. OMBRELO™ offers excellent comfort, stability, and maneuverability. The adjustable, swing-away multi footplates allow easy exit, with foot straps for two stabilization options. Available in 4 sizes, with size 5 supporting adults up to 120 kg.



Lightweight

24<sup>+</sup>

24-month warranty



Stable construction



Easy operation

## Colors

**COMFORT™ upholstery**  
Durable and comfortable



Aloe



Loft



Moon



Onyx



Sand



Rose



Sea



Sky

**HYDRO™ upholstery**  
Waterproof and durable



BodyPrint™ upholstery was equipped with unique belts system for seat and backrest. When adjusted perfectly gives the optimal positioning and comfort of use.



Thanks to multi footplates Flip Flap™ with abduction-adduction angle adjustment stroller grows up with child.



**100 - 180 cm**  
Indicative user height



**45 - 120 kg**  
Max. user weight



For **more details** about the product, please visit our website.

**This is a medical device.**  
**Use it in accordance with the user manual or the label.**

### INDICATIONS

Cerebral palsy (CP) – various types and forms | Brain and spinal cord injuries of various etiologies with limb paresis, severe balance disturbances, and coordination impairments | Demyelinating diseases (multiple sclerosis – MS, amyotrophic lateral sclerosis – ALS, and others) | Meningocele | Muscle diseases with paresis and paralysis (including dystrophy, myasthenia gravis, polymyositis) | Genetic syndromes with limb paresis | Degenerative and metabolic diseases of the nervous system | Degenerative and rheumatic joint diseases | Other diseases with limb paresis, paralysis, and musculoskeletal impairments | After injuries, following traumatic brain injuries or spinal cord injuries with limb paresis and paralysis, particularly with quadriplegia (as part of the rehabilitation process or as permanent orthopedic support)

### CONTRAINDICATIONS

**Observation:** Monitor user during standing. Stop therapy if paleness, sweating, fainting, joint overload, or body deformations occur; consult with a doctor or physiotherapist | **Absolute Contraindication:** Critical condition of the user. | **Relative Contraindications:** Resting hypertension > 200/110 mmHg, pulmonary hypertension, heart rhythm disorders, unstable asthma, unstable metabolic issues, moderate heart valve disease. | Conditions like healing fractures or severe osteoporosis that prevent weight-bearing | Incomplete fracture healing (unless approved by a doctor) | Postoperative joint load restrictions | Pain from dislocation, excessive tissue stretch, or pressure; dislocation without pain is not contraindicated | Significant hip/knee flexion or ankle contractures; ensure no excessive pressure or stretch on affected areas.