

# ProCell II

## Medium risk prevention & treatment of wounds up to grade 2.

ProCell II provides all the essential features an efficient dynamic system needs. The pump offers alternating therapy or constant pressure for a stable low pressure surface. The pressure level of the mattress is monitored and the pump will alert the user if there's an irregularity.

The system includes an overlay mattress to be used on top of an existing foam mattress, but can be equipped with an optimal foam. Underneath the mattress there's room to integrate a foam and make it a mattress replacement.



Quick deflation - CPR  
<18 sec.\*



Foambag for  
Mattress replacement

### Safety Standards

Complies with:  
EC Directive 93/42 Medical Devices  
IEC60601-1: 2005 / EN60601-1:2006

Protection against Electric shock:  
Class II , type BF

Protection of Ingress of Liquids:  
IP21

Mode of Operation:  
Continuous

### Product Specifications

Mattress: Overlay / MR with optional foam

Width: 85cm

Length: 200cm

Height: 13cm (as OL)

Patient weight: 150kg max.

Weight system: 10kg

Pump: 230V / 50Hz

Noise: 26dB

Cycle: 12 min

Design, Quality & Service powered by:

**PRIME**  
Medical Solutions