

ProCell III

Medium risk prevention & treatment of wounds up to grade 3.

The ProCell III system includes a digital pump which offers alternating therapy or constant low pressure with an adjustable weight setting. To assist the caregiver or help the user transfer in or out of bed an Autofirm function can be activated, providing a stable surface when required. For safety this function automatically resets to the previous alternating or static pressure setting after 30 minutes.

Underneath the mattress there's room to integrate a foam and make it a mattress replacement.



Foambag for
Mattress replacement

Safety Standards

Complies with:
Regulation (EU) 2017 / 745
for medical devices

Protection against Electric shock:
Class II , type BF

Protection of Ingress of Liquids:
IP21

Mode of Operation:
Continuous

Product Specifications

Mattress: Overlay / MR with
optional foam

Width: 85 or 90cm

Length: 200 or 210cm

Height: 13cm (as OL)

Patient weight: 150kg max.

Weight system: 10kg

Pump: 230V / 50Hz

Noise: 0 - 26dB

Cycle: 10 / 15 / 20 / 25 min

Design, Quality & Service powered by:

PRIME
Medical Solutions