

Aurora – Riser recliner chair



Meet the new Aurora

Hepro has 35 years of experience in making products that improve independence and quality of life. We focus on innovation and development and work continuously to improve our products.

Now we're launching a new generation of lift chairs. Aurora still has the proven functionalities that its predecessor had, but these features are now more user-friendly. Hepro's new lift chairs, Aurora, have excellent sitting comfort and support for the entire body. The back of the chair adapts well to your body. The way in which it raises you into an upright position and the shape of the armrests make the lift chair very user-friendly. The lift chairs provide comfortable support and make assuming an upright position easy. It is easy to adjust the lift chair and convert it from a vertical lift to a tilt lift.

The chair also has a base plate that does not obstruct the heels when you need to pull your legs well in before getting out of the chair to stand up. All chairs in the new Aurora series come with a swivel function. The swivel release has a better grip and is now easier to manage. In addition to the manual swivel function, there is also an electric swivel option for people with reduced strength in their hands.

Aurora is available in five widths from 40 to 60 cm (the number behind the model name refers to the width of the chair between the armrests measured in cm). In addition, we have a good range of accessories that make the chairs very adaptable. The Aurora lift chair can be adapted for anyone with easy customisation such as finding the right neck pillow or support cushion. The adjustment options allow the chair to be adapted to your needs and condition, e.g. to assume a particular position for eating, resting or sleeping. The Tilt model tilts safely back into a resting position, but also fully back into the Trendelenburg position, i.e. the legs are higher than the head.

Thanks to a USB charger in the chair, phones, tablets or other assistive technology aids are readily accessible, even while charging. Aurora chairs can fit any home and will seem timeless. The materials are made of fabric, leather and oak and are similar to those used in regular furniture. In this way, the chair becomes a normal part of the living room.

We are proud that 30 years later the lifting chairs continue to be manufactured in Rognan in Norway.





All models except 60 cm can be easily converted from a tilt lift to a vertical lift (see the user manual).

Fabric or leather?

The advantage of a fabric lifting chair is that it does not get as slippery and sticky. Furthermore, there is less sliding in a fabric chair. The fabric (excluding the armrests) is removable and can be washed in the washing machine.

It is easier to keep a leather chair clean. It may feel more exclusive to look at, but is a little stiffer and may not feel as soft. We recommend that you use a leather lifting chair if you spill a lot. However, if you suffer from incontinence, we recommend incontinence covers; see pages 8–9.

Aurora

Range of products




Aurora Standard

With individual adjustment of the back and footrest.

Suitable for anyone who needs a lift chair.

The remote control for the Standard chair



User controls back angle/lift and footrest — independently of each other, with 4 buttons.



Aurora Tilt

These models are the most advanced chairs in our range. The chair can be tilted all the way back in a supine position for optimal relief. All of our Tilt models can be tilted back into the Trendelenburg position, i.e. the legs are higher than the head.

The chair is suitable for those who sit for large parts of their day and need varied relief. It is important to change the sitting position to avoid pressure ulcers.

The remote control for the Tilt chair



User controls back angle/lift, footrest and tilt independently of each other with 6 buttons. Automatic lift by pressing the top left button.



Aurora Synchronous

The footrest and back move synchronously to ensure the least possible strain on the hips when changing position. The chair relieves strain on the legs while watching TV or reading.

This chair is particularly suitable for cognitively weak users.

The remote control for the Synchronous chair



The footboard lifts (tilts) to an angle of 23.4 degrees. The back remains stationary.

The back moves back (tilts) to an angle of 58.4 degrees. This opens the hip angle. The footboard retains its position.

The footboard lifts to the maximum position. The back remains stationary.

The back continues backwards. You can now achieve the maximum resting position.

Aurora

Models


Aurora 40
Aurora 45
Aurora 50
Aurora 55
Aurora 60

Leather/Fabric		Leather/Fabric		Leather/Fabric		Leather/Fabric		Leather/Fabric	
Seat width	User weight	Seat width	User weight	Seat width	User weight	Seat width	User weight	Seat width	User weight
40-50 cm	< 165 kg	45-55 cm	< 165 kg	50-60 cm	< 165 kg	55-65 cm	< 165 kg	60-70 cm	<166-260 kg

The number indicates the seat width in cm

Art. No.	Seat width	Description
36100	40 cm	Aurora Synchronous 40 fabric
36101	40 cm	Aurora Synchronous 40 leather
36102	45 cm	Aurora Synchronous 45 fabric
36103	45 cm	Aurora Synchronous 45 leather
36104	50 cm	Aurora Synchronous 50 fabric
36105	50 cm	Aurora Synchronous 50 leather
36106	55 cm	Aurora Synchronous 55 fabric
36107	55 cm	Aurora Synchronous 55 leather
36108	40 cm	Aurora Standard 40 fabric
36109	40 cm	Aurora Standard 40 leather
36110	45 cm	Aurora Standard 45 fabric
36111	45 cm	Aurora Standard 45 leather
36112	50 cm	Aurora Standard 50 fabric
36113	50 cm	Aurora Standard 50 leather
36114	55 cm	Aurora Standard 55 fabric
36115	55 cm	Aurora Standard 55 leather
36116	60 cm	Aurora Standard XXL 60 fabric
36117	60 cm	Aurora Standard XXL 60 leather
36118	40 cm	Aurora Tilt 40 fabric
36119	40 cm	Aurora Tilt 40 leather
36120	45 cm	Aurora Tilt 45 fabric
36121	45 cm	Aurora Tilt 45 leather
36122	50 cm	Aurora Tilt 50 fabric
36123	50 cm	Aurora Tilt 50 leather
36124	55 cm	Aurora Tilt 55 fabric
36125	55 cm	Aurora Tilt 55 leather
36126	60 cm	Aurora Tilt XXL 60 fabric
36127	60 cm	Aurora Tilt XXL 60 leather

Neck pillows



Neck pillow with polystyrene beads

Fabric or leather

Large neck pillow filled with polystyrene beads in 2 compartments. Very soft, comfortable neck pillow that provides good support and shapes itself to the neck and head. The level of softness and support can be changed by altering the amount of polystyrene beads. Polystyrene beads emit heat to tight muscles. Height: 19 cm, Width: 35 cm, Depth: 9 cm.

Art. no. Fabric	Art. no. Leather
36365	36366



Neck pillow with down

Fabric or leather

Large neck pillow with down. Height: 45 cm, Width: 45 cm.

Art. no. Fabric	Art. no. Leather
36373	36374



Neck pillow, dogleg

Fabric or leather

Height: 20/23 cm, Width: 33 cm.

Art. no. Fabric	Art. no. Leather
36304	36303



Neck pillow, medium

Fabric or leather

Foam neck pillow with medium droplet section. Can be used either way to make you comfortable. Height: 18 cm, Width: 36 cm, Depth: 6.5 cm (drop).

Art. no. Fabric	Art. no. Leather
36367	36368



Neck pillows



Extra-large neck pillow (polystyrene beads)

Fabric or leather

Large neck pillow filled with polystyrene beads. Comes with a fabric cover. Very soft and comfortable neck pillow which provides good support. The level of softness and support can be changed by altering the amount of polystyrene beads. Polystyrene beads emit heat to tight muscles. Particularly suitable for people with kyphosis. H: 25 cm, W: 40 cm, D: 12 cm.

Art. no. Fabric	Art. no. Leather
36369	36370



Extra-large (padded) neck pillow

Fabric or leather

Large, soft padded cushion for extra head and neck support. Particularly suitable for people with kyphosis. Height: 25 cm, Width: 40 cm, Depth: 12 cm.

Art. no. Fabric	Art. no. Leather
36371	36372



Stable side-support neck pillow

Fabric or leather

Stable mounted upper back section with soft side supports (ear pads). This provides excellent and stable side support for the head, yet the support is also soft and comfortable. Depth: 12 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36607	36608
45 cm	36609	36610
50 cm	36375	36376
55 cm	36611	36612
60 cm	36613	36614



Accessories

Covers



Super soft back and seat cover

Fabric or leather

A special cover that is pulled over the existing cover with Dacron sewn into the entire cover. Gives even softer sitting comfort. Thickness 3 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36377	36378
45 cm	36379	36380
50 cm	36381	36382
55 cm	36383	36384
60 cm	36385	36386



Super soft back and seat cover + 5 cm

Fabric or leather

A separate cover that is pulled on over the existing cover. The cover is intended for those with deeper seats. It is divided into two, one cover for the back and one for the seat. It is attached with Velcro around the back frame.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36387	36388
45 cm	36389	36390
50 cm	36391	36392
55 cm	36393	36394
60 cm	36395	36396



Backrest comfort

Fabric or leather

Separate cover that is pulled over the lower back section. The cover has 3 large, full-width transverse channels. The channels are filled with Dacron. This provides a soft and comfortable back that is very suitable for people with pain and stability problems.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36407	36408
45 cm	36409	36410
50 cm	36411	36412
55 cm	36413	36414
60 cm	36415	36416



Accessories

Covers



Backrest, side support

Fabric or leather

For extra good side support. A cover that is pulled over the lower back section. The cover has 2 longitudinal wedge-shaped foam channels. This is to provide extra good side support. Two centimetres of Dacron is sewn in between the wedges. Very suitable for people with stability problems. The wedges start 10 cm from each side of the back and are 10 cm deep.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36397	36398
45 cm	36399	36400
50 cm	36401	36402
55 cm	36403	36404
60 cm	36405	36406



Incontinence seat and back cover

Fabric or leather

Incontinence seat and back cover with an incontinence barrier sewn inside. The cover is easy to remove and clean. Comes with attached washing instructions.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36250	36251
45 cm	36252	36253
50 cm	36254	36255
55 cm	36256	36257
60 cm	36258	36259



Incontinence seat cover

Fabric or leather

Incontinence seat cover with an inner incontinence barrier sewn in. The cover is easy to remove and clean. Comes with attached washing instructions.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36260	36261
45 cm	36262	36263
50 cm	36264	36265
55 cm	36266	36267
60 cm	36268	36269



Armrest protection

Fabric or leather

A removable and washable cover that can be put on over your original armrests for increased protection. This makes cleaning easier and improves hygiene.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36246	36247

Accessories

Seat cushions



Seat cushion with pressure relief

Fabric or leather

With pressure-relieving foam.

Depth: 52 cm, Thickness: 11 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36447	36448
45 cm	36449	36450
50 cm	36451	36452
55 cm	36453	36454
60 cm	36455	36456

Gel seat cushion

Fabric or leather

Is made of foam and gel pockets. A softer cushion for increased pressure relief to avoid pressure ulcers.

Depth 52 cm, Thickness 11 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36520	36521
45 cm	36522	36523
50 cm	36524	36525
55 cm	36526	36527
60 cm	36528	36529



Deep seat cushion

Fabric or leather

A deeper seat cushion for those who find the standard cushion too short. Depth 55 cm, Thickness 11 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36437	36438
45 cm	36439	36440
50 cm	36441	36442
55 cm	36443	36444
60 cm	36445	36446



Low seat cushion

Fabric or leather

Half-thickness seat cushion (5 cm) for a lower sitting height.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36427	36428
45 cm	35429	36430
50 cm	36431	36432
55 cm	36433	36434
60 cm	36535	36436

Accessories

Seat cushions



Small flat seat

Leather

The flat seat can be used either way, so that the edge is in front or behind.

Depth: 45 cm.

Width	Art. no.
40 cm	36457
45 cm	36458
50 cm	36459
55 cm	36460
60 cm	36461



Large flat seat

Leather

The flat seat can be used either way, so that the edge is in front or behind.

Depth: 50 cm.

Width	Art. no.
40 cm	36462
45 cm	36463
50 cm	36464
55 cm	36465
60 cm	36466



Accessories

Back cushions/support



Lumbar support, small (inflatable)

Fabric or leather

Inflatable air cushion with a leather or fabric cover and a separate hand pump. The cushion can be attached anywhere to the back of chair as needed with a belt. 20 x 40 cm.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36423	36424



Lumbar support, large (inflatable)

Fabric or leather

Inflatable air cushion with a leather or fabric cover and a separate hand pump. The cushion can be attached anywhere to the back of chair as needed with a belt. 25 x 50 cm.

Art. no. Fabric	Art. no. Leather
36425	36426



Back cushion, small

Fabric or leather

Foam cushion with a triangular cut. The cover is shaped so that the cushion can be secured around the back of the chair. The cushion can be attached anywhere to the back of chair as needed with a belt. Height: 19 cm.

Width	Art. no. Fabric	Art. no. Leather
40-45 cm	36417	36418
50-55 cm	36419	36420
60 cm	36421	36422



Accessories

Footrest



Footrest, extra wide

Fabric or leather

Width: 70 cm.

Art. no. Fabric	Art. no. Leather
36477	36478



Footrest with side support

Fabric or leather

4.5 cm high side supports.

Width	Art. no. Fabric	Art. no. Leather
40 cm	36479	36480
45 cm	36481	36482
50 cm	36483	36484
55 cm	36485	36486
60 cm	36487	36488



Accessories

Belts



Belt restraint

Small

Child size.

Art. no.
27850



Belt restraint

Medium

Smaller adult size.

Art. no.
27851



Belt restraint

Large

Adult size.

Art. no.
27852



Belt restraint

Exstra Large

Large adult size.

Art. no.
27853



Hip belt

Small

Hip belt that provides effective and comfortable support and improves pelvic stability.

Art. no.
27847



Hip belt

Medium

Hip belt that provides effective and comfortable support and improves pelvic stability.

Art. no.
27848



Hip belt

Large

Hip belt that provides effective and comfortable support and improves pelvic stability.

Art. no.
27849



Magazine pocket

Fabric or leather

Spacious pocket that hangs over the armrest on the right or left side of the chair. For storing newspapers, magazines and other things you want to keep close by.

Art. no. Fabric
36248

Art. no. Leather
36249



Crutch holder

Can be mounted on both sides.

Art. no.
26565

Accessories

Miscellaneous



Round table

Small round table for, maybe, a coffee cup. Fits under the armrest on the right or left side. Diameter 27.5 cm.

Art. no.
26591



Rectangular table

31x41 cm

A small table that can be fitted to the right or left side of the chair. Mounted on the frame of the chair. Perfect if you want to read, write or use a laptop in the chair.

Side	Art. no.
Left	27839
Right	26592



Extension set, foot lift scissors

A separate scissor joint that is easily mounted into the foot-lift mechanism so that the footrest comes out a further 13 cm. Includes fabric or leather joints.

Art. no.
36494
36495



Handset

Modell	Art. no.
Synchronous (2 buttons)	61817
Standard (4 buttons)	61816
Tilt (6 buttons)	61815



Button raiser

Transparent rubber button that is placed over the remote control. This raises the buttons and makes them easier for the visually impaired to operate. Comes in a packet with six pieces.

Art. no.
29640



Joystick control

The joystick is mounted on its own adjustable bracket next to the armrest on the right or left side. This can control 2 motors. For a 3-motor chair, the last motor must be controlled using a separate switch.

Art. no.
28684 (standard)
28909 (tilt)

Accessories

Miscellaneous



Battery backup

If there is no access to electricity, e.g. at the holiday cabin or with an unstable power supply, a battery backup can be used. This has rechargeable batteries that allow the chair to rise and lower a few times so that it does not get stuck in a particular position.

Art. no.

28683



Ball in front of wooden armrest

A ball in the front of the armrest for increased grip when sitting down on and getting up from the chair.

Art. no.

80007



Castor chassis

With central braking

Castor chassis for your Aurora lift chair which makes it easier for users to move around without having to leave the chair. The chassis is supplied with a brake lever.

Art. no.

Fits

27723 Castor chassis with central braking Aurora XXL Tilt

26524 Castor chassis with central braking Aurora XXL

29641 Castor chassis with central braking Aurora Tilt

29085 Castor chassis with central braking Aurora



Electrical rotation

Electrical rotation of the chair by remote control. Suitable for people with reduced strength or functional ability of their legs.

Width

Art. no.

Fits

40-55 cm

29993

Standard

40-55 cm

29984

Tilt

60 cm

29903

Standard

60 cm

29901

Tilt



Measurements and dimensions

	Model 40	Model 45	Model 50	Model 55	Model 60
Front seat height min.	45 cm	45 cm	45 cm	45 cm	47 cm
Front seat height max.	55 cm	55 cm	55 cm	55 cm	57 cm
Seat height in upright position, tilt lift, front	59-69 cm	59-69 cm	59-69 cm	59-69 cm	59-69 cm
Seat height in upright position, tilt lift, rear	71-78 cm	71-78 cm	71-78 cm	71-78 cm	73-80 cm
Seat height in upright position, vertical lift, front	69-79 cm	69-79 cm	69-79 cm	69-79 cm	-
Seat height in upright position, vertical lift, rear	73-80 cm	73-80 cm	73-80 cm	73-80 cm	-
Seat depth min.	36 cm	39 cm	39 cm	39 cm	42 cm
Seat depth max.	45 cm	48 cm	48 cm	48 cm	52 cm
Seat width	40 cm	45 cm	50 cm	55 cm	60 cm
Back height	81-99 cm	81-99 cm	81-99 cm	81-99 cm	81-99 cm
Armrest distance min. (inner)	40 cm	45 cm	50 cm	55 cm	60 cm
Armrest distance max.	50 cm	55 cm	60 cm	65 cm	70 cm
Armrest height min.	6 cm	6 cm	6 cm	6 cm	6 cm
Armrest height max.	22 cm	22 cm	22 cm	22 cm	22 cm
Armrest width	13 cm	13 cm	13 cm	13 cm	13 cm
Largest width of chair	67 cm	72 cm	77 cm	82 cm	87 cm
Maximum length of chair	182 cm	182 cm	182 cm	182 cm	182 cm
Tilt in degrees (for tilt chairs)	23 °	23°	23°	23°	23°
Max. user weight	165 kg	165 kg	165 kg	165 kg	260 kg
Chair weight	57 kg	59 kg	61 kg	63 kg	79 kg

Aurora

Standard parts

Neck pillow

Fabric

Width	Art. no.	Dimension
40 cm	36467	H: 23 cm, W: 30 cm
45 cm	36468	H: 23, cm, W: 35 cm
50 cm	36281	H: 23 cm, W: 40 cm
55 cm	36469	H: 23 cm, W: 45 cm
60 cm	36470	H: 23 cm, W: 50 cm

Neck pillow

Leather

Width	Art. no.	Dimension
40 cm	36238	H: 23 cm, W: 30 cm
45 cm	36239	H: 23, cm, W: 35 cm
50 cm	36336	H: 23 cm, W: 40 cm
55 cm	36240	H: 23 cm, W: 45 cm
60 cm	36241	H: 23 cm, W: 50 cm

Seat cushion

Fabric

Width	Art. no.	Dimension
40 cm	36274	Depth: 50 cm
45 cm	36286	Depth: 52 cm
50 cm	36296	Depth: 52 cm
55 cm	36308	Depth: 52 cm
60 cm	36318	Depth: 52 cm

Seat cushion

Leather

Width	Art. no.	Dimension
40 cm	36332	Depth 50 cm
45 cm	36340	Depth: 52 cm
50 cm	36347	Depth: 52 cm
55 cm	36353	Depth: 52 cm
60 cm	36359	Depth: 52 cm

Aurora

Standard parts

Footrest cushion

Fabric

Width	Art. no.	Dimension
40 cm	36283	Height: 27 cm
45 cm	36293	Height: 27 cm
50 cm	36305	Height: 27 cm
55 cm	36315	Height: 27 cm
60 cm	36325	Height: 27 cm

Footrest cushion

Leather

Width	Art. no.	Dimension
40 cm	36338	Height: 27 cm
45 cm	36344	Height: 27 cm
50 cm	36351	Height: 27 cm
55 cm	36357	Height: 27 cm
60 cm	36363	Height: 27 cm

Additional customisation of your Aurora chair

Neck problems	- Adjust the angle of the upper back section to the desired position - Choose a neck pillow of the right size and softness. See examples in the list of accessories
Balance problems and need for positioning in the chair	Choose one of our back covers for extra support. For example, a back cover with extra lateral support
Need for extra adjustment of the back, e.g. for kyphosis, scoliosis and similar issues	Use, for example, a Comfort back cover with fibre filling that can be shaped according to the user's needs
Impaired function of the fingers and/or hands	Use the joystick instead of the handset.l.
Need for extra soft sitting comfort	- Use the Supersoft cover - Use the Comfort back cover - Use the pressure-relieving seat cushion
Need for extra lumbar support	- Use the inflatable lumbar support. The cushion is available in two different sizes. - If necessary, use the small back cushion, which can be attached to the back as required
Incontinence problems	Use the incontinence cover for the seat, seat/back and, if applicable, the incontinence cover for the footrest.
Risk of pressure ulcers	Use the pressure-relieving seat cushion, the gel seat cushion, the Supersoft cover or a flat seat when using special cushions
Need for support for longer legs	Use the scissor mechanism lift extension set. The footrest will then come forward a further 13 cm
Need to operate chair from the left	The swivel lever and the handset holder can be moved to the left side on request. Can also be retrofitted.
Risk of falling out of the chair	Use belt holder and/or hip belt
Need to move the chair without getting out of it	Use the castor chassis with central braking

Errors and omissions excepted.



Active and safe at every stage of life

HEPRO AS is a modern and future-oriented competence company that develops, produces and sells technical aids for children, young people and adults, along with high quality assistance technology. Our head office is located in Rognan, at the heart of Saltenfjorden, surrounded by mountains and magnificent scenery.

The company was established in 1987 and has extensive experience in developing products and solutions to make everyday life easier for our customers. Our vision is to enable our customers to manage on their own, actively participate in their own everyday life and break down barriers.

Through high quality, fast response times and a short communication chain, Hepro shall be the natural first choice for our customers.

In close cooperation with occupational therapists, we develop products and solutions tailored to our customers' needs for help with mobility, balance, comfort and hygiene. Our aim is to provide confidence, the ability to cope and increased activity through functional and user-oriented products and technology. We engage with the individual user and go above and beyond to deliver products tailored to their individual needs.