

Roger™ receiver overview

(Options 02 and 03)

Receivers colors

roger

Color overview



Roger 10
AS10



Roger 11
AS11



Roger 14



Roger 15
AS15



Roger 16
AS16



Roger 17



Roger 18
AS18



Roger 19
AS19



Roger 21

		Roger 10 AS10	Roger 11 AS11	Roger 14	Roger 15 AS15	Roger 16 AS16	Roger 17	Roger 18 AS18	Roger 19 AS19	Roger 21
Sand Beige	P1									
Amber Beige	P2									
Sandalwood	P3									
Chestnut	P4									
Champagne	P5									
Silver Gray	P6									
Graphite Gray	P7									
Velvet Black	P8									
Ruby	P9									
Petrol	Q1									
Electric Green	Q2									
Caribbean Pirate	Q3									
Dragon Orange	Q4									
Vanity Pink	Q5									
Blue Lagoon Transparent	Q6									
Royal Purple Transparent	Q7									
Precious Pink	T3									
Beige	01									
Brown	03									
Light Gray	04									
Black	06									
Pure Transparent	13									
Blue Transparent	37									
Purple Transparent	38									
Gray	23									
Light Pink	33									
Light Blue	34									
Taupe	35									
Chestnut Brown	50									
Palladium Black	70									
Chocolate	E5									
Cliff	E6									
Mid Gray	E8									
Glacier	E9									
Dark Cliff Gray	F0									
Sterling	F3									
Brown	L0									
White	L8									
Black	L9									
Beige	M1									
Charcoal	M2									
Mocha	T1									
Smoke	T2									
Alpine White	XN/T7									
Princess Pink	XP									
Lava Red	M6									
Blue Ocean	M7									
Majesty Purple	M8									
Anthracite	V1									
Beige	V2									
Black	V3									
Ebony	V4									
Nordic Grey	V5									
White	V6									

Roger receiver compatibility guide – Hearing aids

		Design-integrated for Phonak hearing aids						Universal			
		Roger 19	Roger 18	Roger 16	Roger 15	Roger 11	Roger 10	Roger X	Audio shoe	Roger MyLink	
Belong	RIC	Audéo B-13		•					•	AS18	•
		Audéo B-312T							• ¹		•
		Audéo B-312							• ¹		
		Audéo B-10							• ¹		
		Audéo B-R							• ¹		
		Bolero B-M							• ¹		•
	BTE	Bolero B-P		•					•	AS18	•
		Bolero B-SP		•					•	AS18	•
		Bolero B-PR							• ¹		
		Phonak Sky V-M							• ¹		•
	Pediatric	Phonak Sky V-P		•					•	AS18	•
		Phonak Sky V-SP		•					•	AS18	•
Phonak Sky V-UP		•						•	AS19	•	
Phonak Sky V-RIC			•					•	AS18	•	
Phonak Naida V-SP			•					•	AS18	•	
Phonak Naida V-UP		•						•	AS19	•	
Venture	Power	Phonak Naida V-RIC		•				•	AS18	•	
		Phonak Bolero V-M							• ¹		•
		Phonak Bolero V-P		•					•	AS18	•
		Phonak Bolero V-SP		•					•	AS18	•
	RIC	Audéo V-13		•					•	AS18	•
		Audéo V-312T							• ¹		•
		Audéo V-312							• ¹		
		Audéo V-10							• ¹		
		Virto V-10 O							• ¹		•
		Virto V-10							• ¹		•
	Custom	Virto V-312							• ¹		•
		Virto V-13							• ¹		•
Phonak Sky Q-M13					•			•	AS15	•	
Phonak Sky Q-SP						•		•	AS11	•	
Quest	Pediatric	Phonak Sky Q-UP					•	•	AS10	•	
		Phonak Sky Q-RIC				•			•	AS15	•
		Phonak Naida Q-SP					•		•	AS11	•
		Phonak Naida Q-UP						•	•	AS10	•
	Power	Phonak Naida Q-RIC				•			•	AS15	•
		Phonak Bolero Q-M 312							• ¹		•
		Phonak Bolero Q-M 13				•			•	AS15	•
		Phonak Bolero Q-P			•				•	AS16	•
		Phonak Bolero Q-SP			•				•	AS16	•
		Phonak Baseo Q-M 13				•			•	AS15	•
	Basic	Phonak Baseo Q-P			•				•	AS16	•
		Phonak Baseo Q-SP			•				•	AS16	•
Audéo Q 312T								• ¹		•	
Audéo Q 312								• ¹			
Custom	Virto Q 10										
	Virto Q 312							• ¹		•	
	Virto Q 13							• ¹		•	

1 Roger X must be used with a ComPilot or ComPilot II streamer

Some older Spice H20 models are compatible with Roger 15

Roger receiver compatibility guide – Implants

		Design-integrated for implants			Universal			
		Roger 14	Roger 17	Roger 21	Roger X	Audio shoe	Roger MyLink	
Implants	AB	Naida CI Q		•		• ¹	•	
		Harmony / Auria				•	iConnect	•
		Neptune				•	Neptune Connect	•
	Cochlear	Nucleus 5	•			•	Euro accessory adapter	•
		Nucleus 6	•			•	Euro accessory adapter	•
		Baha 5				• ²		
		Baha 4				•		•
		Baha BP100 / BP110				•		•
		Freedom						•
		Sonnet			•	•	FM battery pack cover	•
	Med EI	Rondo				•	Mini battery pack	•
		Opus2				•	FM battery pack cover	•
	Oticon Medical	Ponto Plus / Plus Power				•	Streamer	•
		Ponto / Ponto Pro / Ponto Pro Power						• ³

1 Roger X must be used with a ComPilot or ComPilot II streamer

2 Via audio cable connected to Cochlear Mini Microphone

3 Requires external Oticon medical Telecoil

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we remain true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonakpro.com