

Invacare® | *Alu Lite*



A lighter way to get around

The *Invacare Alu Lite* is the basic-entry wheelchair in the Invacare manual range. The lightweight aluminium construction of this transit wheelchair makes it easy to manoeuvre, even for an elderly carer, as well as providing a wheelchair that is easily transportable.



Yes, you can.®



Features and Options

Easy to transport and store

Being lightweight (aluminium chassis), compact and with a foldable backrest, the **Invacare Alu Lite** is very easy to handle, transport and store. Simply remove the footrests, fold the backrest and chassis and store away.

Pushing? No problem!

The overall manoeuvrability of a wheelchair is crucial, which is why the **Alu Lite** is designed to be very easy to push and manoeuvre, both indoor and outdoor. Due to the smaller overall dimensions and general compactness of this model, getting around in tight spaces is never a problem.



Folding backrest

Reduces overall size when folded for transportation or storage.



Swing away footrest









Footplates fold up as shown and footrests swing away, allowing easier access in or out of the chair.














Foldable cross frame

Lightweight frame folds flat for easier portability.

Technical data

								
Invacare Alu Lite	405 / 455 mm	430 mm	0°	490 mm	0°	450 mm	80°	340 - 460 mm

						
Invacare Alu Lite	230 mm	Seat width + 110 mm	250 mm	1030 mm	760 mm	945 mm

					
Invacare Alu Lite	12,5°	890 mm	100 kg	12 kg	10,2 kg

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 GL
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeurope@invacare.com
www.invacare.eu.com

All product dimensions are in millimetres to apply EN 12183 standards. For more pre-sales information about this wheelchair, please consult the www.invacare.eu.com website. It can be found under the section Low Active Wheelchairs / Alu Lite / Sales and Images in the product catalogue. For information about typical user description, intended use, adjustments and transportation, please refer to the user guide.

Invacare International Sàrl
Route de Cité-Ouest 2 - 1196 Gland
Switzerland
Tel: +41 (0)22 354 60 10
Fax: +41 (0)22 354 60 11
hqeuropa@invacare.com
www.invacare.eu.com

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Alu Lite - EU - 06/2011

Frame colour



Blue

Upholstery



Black nylon



Please notice colors may vary slightly from those displayed above.



Yes, you can.®