CLIP







Functions





Net basket on the back of the backrest (pict.1)

The CLIP is **light weighted** and can be **easily folded**

The CLIP has an **inclined seat** to assure the **maximum comfort** to the child.

Folding aluminium frame







Accessories

812 Directional locks



▶ 827 Foot straps





▶ 816 Side supports for the trunk



▶ 817 Leg straps for abduction



▶ 835 adjustable Headrest



894 Pelvic belt



853 Vest harness



Canopy



▲ 818 Thermal cover



Five point harness

906



819



▲ 825 Rain cover to be fixed to the 819 canopy



WARNING:


- fits the ISO 10542 compliant 4-POINT STRAP-TYPE TIEDOWN, provided by the vehicle
- it is required for the child to wear both pelvic- and shoulder-belt-type restraints, to be compliant either with ISO 7176-19 or ISO 10542
- postural belts and harnesses such as our 828, 853, 894, 903, 906 are not approved for the use as an occupant restraint system
- the vehicle must be suitable for transporting orthopedic pushchairs in compliance with current law
- read the warnings and instructions given in the user's & maintenance book



Five point vest harness



891 SET OF 4 TIE DOWN HOOKS (only SIZES 1-2-3) for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses..)



THE NUMBERS DIVIDED BY THE DASH SPECIFY A MINIMUM AND A MAXIMUM ADJUSTMENT

Product Certified by the following institute: Product marked Product

 ϵ

FDA USA

registered

Authorized Dealer

CERAH France

