



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE ALMOST FOUR YEARS I'VE HAD IT. IT HAS NEVER BROKEN DOWN ON ME OR LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays - Texas, USA

magic
Mobility™

Frontier V6 All Terrain
**FOR ULTIMATE
OUTDOOR
ADVENTURE**

FULL LIGHTING PACKAGES

Been out having too much fun and it's getting late? Put your lights on. Fitted on rubber brackets.



UP CLOSE AND PERSONAL

Side transfers aren't the only reason you need to flip your armrest out of the way.



REHAB SEATING

Infinite options for infinite body shapes. Our rehab seating range is easily adjustable or can be custom designed to a fixed position. It's whatever you need.

TAKE CONTROL

A full range of advanced control systems and environmental controls are available. Easy to program and easy to use.



Frontier
V6™

V3/AUST/MARCH19

BUILT TO LAST

CUSTOMISED FOR YOUR LIFESTYLE



AUSTRALIAN MADE



FRONTIER V6 ALL TERRAIN

Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer the track to the tarmac.

Basic specifications:**

- Length: 1050mm
- Width: 710mm
- Seat to floor height: 435mm
- Turning radius: 580mm
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and tilt
- Top speed: 10kph (6.2mph)
- Motors: 800W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Castors: 210x60mm (2.50-4) black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-24" wide)
- WC19 crash test compliant tie down points.

Rehab positioning:

- Power seat options: 12" lift, tilt (45 degrees), lift and tilt, transfer tilt with tilt only, transfer tilt with lift and tilt, power recline, lay down, power elevating legrests
- Seating and Backrests: MPS, rehab backrest, adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: Centre post, swingaway, power and manual elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Retractable docking pin
- Puncture resistant tyre liners (off-road wheels only)
- 24V LED Lighting kit
- USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, drinks holder, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
- Hybrid wheel kit: narrower, grey and black drive wheels and castors, when off-road capabilities are not needed. This also reduces the size of the chair further
- RAM mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

The right wheelchair for:*

- Adventurous off-road recreation
- Everyday off-road, country living and working
- Urban outdoor living and working
- Roomy indoor spaces.

Base options:



V6 ALL TERRAIN - BLACK RIMS



V6 HYBRID - BLACK WHEELS



V6 HYBRID - GREY WHEELS

"WE'VE BEEN
CREATING CUSTOM
MADE VEHICLES
FOR 25 YEARS,
SO THAT YOU CAN
LIVE THE LIFE
YOU CHOOSE"

Ashley Daff
Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.

BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.



Magic Mobility +61 3 8791 5600
sales@magicmobility.com.au

Explore with us magicmobility.com.au



Scan for agent locations

*Please remember that your wheelchair has electric motors and must NEVER be driven through water, into rivers, creeks and the sea

**Specifications subject to change. Check online for current version

**Various model options excluded