



**Gel
Ovation**

Pressure Care Solutions

PRODUCT SELECTOR



Life changing

Gel Oventions products help maintain skin integrity, improve posture and reduce the risk of pressure ulcers.



About us

Established in 2007, Gel Oventions products are used globally throughout healthcare sectors, including the NHS, providing pressure management solutions for wheelchairs, mobility aids and all types of seating. Widely issued by Occupational Therapists,

Physiotherapists, Rehabilitation Specialists and Nurses, our family of Gel-based medical products provide over 300 solutions for your pressure care requirements.

Our Dimensional Gel Pads

Our Dimensional Gel Pads are an essential part of pressure care management. Used and trusted throughout the world in hospitals, care homes and organisations including the NHS and Community Stores. The pads provide a soft, cooling and comforting layer of protection for the prevention of pressure ulcers.



Pressure Protection

Our unique Dimensional Gel Pads help maintain skin integrity and reduce the risk of pressure ulcers.

Wheelchair Accessories

Our gel based wheelchair accessories are engineered to ensure the ultimate in comfort, flexibility and locked in strength.



Bathroom Accessories

Gel Oventions Shower and Commode Accessories help protect the skin when it is at its most vulnerable. These easy to fit options help contribute to a 24 hour skin protection plan.

Products and their Benefits

Soft Cool Gel

All Gel Oventions products incorporate a soft cool gel in various forms that provide tried and tested pressure solutions. The products are designed to protect the skin against surfaces that are hard and unyielding and offer little or no protection to a user's skin. We offer comfort and protection whenever and wherever contact is unavoidable.

Gel Wheelchair Accessories and Hardware

Instant Solutions

When pressure is an issue, people often need to find an instant solution. That is why all our products are designed to be Easy Fit. Most are either Stick On, Velcro On, or can be Clamped or Bolted into place with simple brackets or bolts.

Ergonomic Design

All our products are ergonomically designed to look and feel great. We want them to appear as if they are an original part of the wheelchair or mobility aid, but also add extra comfort that is better than the original. Gel Oventions products avoid unsightly adaptations on otherwise fantastic looking mobility aids.

Secure And Sturdy Products

Our high performance hardware accessories are constructed from high-strength aluminium alloy and are engineered for strength and flexibility.

Our hardware supports feature powerful locking mechanisms for which Gel Oventions has become renowned. This keeps the unit securely in place for maximum user safety and stability at all times.

Economical

By providing an instant solution, our products can eliminate the need for a repeat clinical assessment or home visit. This saves considerable costs in time and transport whilst also immediately helping the end user by eliminating the pressure risk at an early stage.

Many modifications can take long hours of labour and expensive materials to create a suitable solution whilst the user is left waiting. We offer products that can go straight from the shelf to the user saving valuable time, money and discomfort.



Easy Maintenance

Our Dimensional Gel Pads can be easily cleaned using a 10:1 water to bleach solution to sanitize the pad. For NHS establishments please use your Trust approved cleaning agent, then rinse with water and dry with a towel. The pads can also be supplied with a durable VP Cover.



Throughout this catalogue, products with this symbol contain our Dimensional Gel.



The Adjustable Forehead Stabiliser

See page 14



This multi-positional unit is designed to improve the head and neck position enhancing posture, vision, communication skills and daily living activities.

Lateral Supports and Hardware

See page 42



Gel Oventions make a variety of products designed to offer the client comfort and protection whenever and wherever contact is inevitable.

Arm and Elbow Supports

See page 30




We have many easy fit options for Arms and Elbows that provide positioning combined with the comfort of gel.



Leg and Foot Supports

See page 19



From simple Gel Wraps to Multi Axis Hip Guides Gel Oventions have numerous options to support lower limbs.

Wraps and Covers

See page 47

Dimensional Gel Pads

Gel Ovation's Dimensional Pads are made of medical-grade Silicone Polymer designed specifically to have the same physical characteristics as human fat tissue. The softness of the gel stays unchanged at any temperature and when cut, the gel will not flow or leak. It will not dry out like foam. The cellular design allows for even pressure distribution and contours to any underlying surface.

After

The cellular design allows air to circulate between the cells

- minimises heat and moisture build-up
- allows for even pressure
- contours to underlying surface
- drastically reduces shear and friction forces on the skin



Before



correct seat to floor and armrest height

"Rather than replacing Terry's armchair to incorporate pressure care, I prescribed a Gel Ovation's Dimensional Gel Pad to meet his pressure care needs and maintain good posture. This was a cost effective and timely solution."

Danielle - Specialist Occupational Therapist



KEY BENEFITS:



- Prevents Pressure Ulcers
- Expedites Ulcer Healing
- Low Profile
- Improves comfort
- Reduces shear and abrasions
- Improved cooling
- Creates Air-flow
- Hypoallergenic
- Easy to maintain
- Odourless
- Unique Cellular Design
- Contours to any Shape
- Temperature Resistant
- Easy to Sanitise

TYPICAL USES:



- Wheelchair Seating
- Bespoke Seating Systems
- Rise and Recline Chairs
- High Back Chairs
- Beds
- Manufacturing



PRESSURE MAPPING



Client A- Sat in Arm Chair with no pressure protection



Client A- Sat in Arm Chair with the addition of the Dimensional Gel Pad

Good results have been achieved from Pressure Mapping the Dimension Gel Pads. However, the Gel Pads have additional benefits that cannot be recorded through traditional pressure mapping. These include the anti-shear and friction properties, airflow and coolness of the Gel.

We manufacture many comfort and positioning products using our Dimensional Gel, please visit www.gelovationseurope.com

LARGE DIMENSIONAL GEL PADS 1/2" (13MM) THICK

14145D	14" x 14" (360mm x 360mm)
16165D	16" x 16" (400mm x 400mm)
16185D	16" x 18" (400mm x 460mm)
16205D	16" x 20" (400mm x 500mm)
18165D	18" x 16" (460mm x 400mm)
18185D	18" x 18" (460mm x 460mm)
18205D	18" x 20" (460mm x 500mm)
20185D	20" x 18" (500mm x 460mm)
20205D	20" x 20" (500mm x 500mm)
22185D	22" x 18" (560mm x 460mm)
22205D	22" x 20" (560mm x 500mm)
24205D	24" x 20" (610mm x 500mm)
26205D	26" x 20" (660mm x 500mm)
28205D	28" x 20" (711mm x 500mm)
30245D	30" x 24" (762mm x 610mm)



DARTEX COVERS FOR DIMENSIONAL GEL PADS

For the ultimate protection of your Dimensional Gel Pad we also provide Dartex Covers that are available in all the above sizes. The vapour permeable fabric is easy to wipe clean, has a non-slip base and allows the Gel Pad to be removed for cleaning.

To order please use the code: LDCNS followed by the size of pad it is intended to be used with. For example to order a cover for an 18 x 18" pad the order code is: LDCNS1818.

"I regularly use the Dartex non-slip covers for the Dimensional Gel Pads to keep them in place on the chair, wheelchair or shower commode".

Danielle - Specialist Occupational Therapist



SMALL DIMENSIONAL GEL PADS ½" (13MM) THICK

445D	4" x 4" (101mm x 101mm)
355D	3" x 5" (76mm x 127mm)
465D	4" x 6" (101mm x 152mm)
45755D	4.5" x 7.5" (114mm x 190mm)
5755D	5" x 7" (127mm x 178mm)
6755D	6" x 7.5" (152mm x 190mm)
895D	8" x 9" (203mm x 229mm)
555D	7" x 7" (127mm x 127mm)
775D	7" x 7" (178mm x 178mm)
65955D	6.5" x 9.5" (165mm x 241mm)
1285D	12" x 8.5" (305mm x 216mm)

Hexagonal Shaped Pad



SPINAL DIMENSIONAL GEL PADS ½" (13MM) THICK

3105D	3" x 10" (76mm x 260mm)
3125D	3" x 12" (76mm x 4300mm)
3145D	3" x 14" (76mm x 360mmn)
3165D	3" x 16" (76mm x 400mm)
3185D	3" x 18" (76mm x 460mm)
3205D	3" x 20" (76mm x 500mm)



"I often provide a Dimensional Gel Pad to clients that have a tendency of sliding or moving in their chair. The Pad reduces the risk of shear and friction and helps to maintain their position and skin integrity".

Danielle - Specialist Occupational Therapist



Flat Gel Pads

Our Dimensional Gel is also available in a flat style that comes in various thicknesses and has many uses. It reduces point pressure, gradually equalizes to body temperature and is made from Medical Grade Silicone. The Gel will not break down and can be easily cut to any shape or size.



FLAT GEL PADS SPECIFICATIONS

SIZE	1/8" (3.175mm)	1/4" (6.35mm)	3/8" (10mm)	1/2" (13mm)	5/8" (15.87mm)
4" x 4" (100mm x 100mm)	18SFP0404	14SFP0404	38SFP0404	12SFP0404	58SFP0404
5" x 6" (130mm x 150mm)	18SFP0506	14SFP0506	38SFP0506	12SFP0506	58SFP0506
7" x 7" (180mm x 180mm)	18SFP0707	14SFP0707	38SFP0707	12SFP0707	58SFP0707
8" x 8" (200mm x 200mm)	18SFP0808	14SFP0808	38SFP0808	12SFP0808	58SFP0808
8" x 10" (200mm x 260mm)	18SFP0810	14SFP0810	38SFP0810	12SFP0810	58SFP0810
10" x 10" (260mm x 260mm)	18SFP1010	14SFP1010	38SFP1010	12SFP1010	58SFP1010
12" x 12" (300mm x 300mm)	18SFP1212	14SFP1212	38SFP1212	12SFP1212	58SFP1212
14" x 14" (360mm x 360mm)	18SFP1414	14SFP1414	38SFP1414	12SFP1414	58SFP1414
16" x 16" (400mm x 400mm)	18SFP1616	14SFP1616	38SFP1616	12SFP1616	58SFP1616
18" x 16" (460mm x 400mm)	18SFP1816	14SFP1816	38SFP1816	12SFP1816	58SFP1816
18" x 18" (460mm x 460mm)	18SFP1818	14SFP1818	38SFP1818	12SFP1818	58SFP1818
20" x 18" (500mm x 460mm)	18SFP2018	14SFP2018	38SFP2018	12SFP2018	58SFP2018
20" x 20" (500mm x 500mm)	18SFP2020	14SFP2020	38SFP2020	12SFP2020	58SFP2020
22" x 20" (560mm x 500mm)	18SFP2220	14SFP2220	38SFP2220	12SFP2220	58SFP2220
24" x 20" (610mm x 500mm)	18SFP2420	14SFP2420	38FSP2420	12SFP2420	58SFP2420

Bespoke Gel Solutions

As the Dimensional Gel can be cut to any shape or size we can also offer bespoke solutions for you or your client's needs. Popular requests include toilet / shower seat padding, foot and leg rest padding and padding for exposed tubes or bolts. This is a great solution when one of our regular products does not fit your requirements.



“Hi I’m Anna from the Customer Services Team. We understand how important it is that our Pressure Care Products reach our customers as quickly as possible. We aim to dispatch all regular orders within 48 Hours. In fact, over 90% of our orders are delivered the next working day!”



The Adjustable
Forehead Stabiliser
– The latest addition
to Gel Ovations
Hardware Solutions



Adjustable Forehead Stabiliser

This multi-positional unit is designed to improve the head and neck position enhancing posture, vision, communication and daily living activities. Constructed from high-strength light-weight aluminium alloy, the Forehead Stabiliser features Dimensional Gel for supreme comfort, coolness and pressure care. The powerful double locking mechanism keeps the unit securely in place for maximum user safety and stability at all times.

Key features

- Adjustability of pad position places it right where needed
- Works with virtually all headrest systems
- Easily flips back 180 degrees when no longer needed
- Pad allows head movements without shearing or pressure
- Stable, exclusive Gel Ovation design
- Holds its position regardless of client activity
- Affordable and easy to clean
- Incorporates Dimensional Gel
- Maximum comfort & skin protection



HOW TO ORDER

The Multi-Adjustable Forehead Stabiliser can be ordered with one of three sizes of Gel Ovations head rests that incorporate Dimensional Gel. All Gel Ovations Head Rests are compatible with a 1" (25mm) Ball fitting.

Additionally, when a more supportive headrest is required, Adaptor Brackets can be ordered to fit the Forehead Stabiliser to 3rd Part Head Rests including Matrix, Jay, Otto Bock, Stealth, Whitmeyer and all major brands. Please see order codes below.

SRAFS510	SRAFS512
Adjustable Forehead Stabiliser	Adjustable Forehead Stabiliser
Supplied with 5 x 10" (130 x 260mm) Head Rest Pad	Supplied with 5 x 12" (130 x 300mm) Head Rest Pad
SRAFI4	SRAFS
Adjustable Forehead Stabiliser	Adjustable Forehead Stabiliser
Supplied with 5 x 14" (130 x 360mm) Head Rest Pad	Hardware + Forehead Pad ONLY To be fitted to 3 rd Party Head Rests



extension available



FOREHEAD STABILISER MOUNTING BRACKETS FOR 3RD PARTY HEAD RESTS

SRFSMB	SRFSMBS
Stabiliser Mounting Bracket and Screws for 1" (25mm) Ball Fitting	Stabiliser Mounting Bracket and Screws for ¾" (19mm) Ball Fitting
SRFSMBW	SRFSMBC4
Stabiliser Mounting Bracket and Screws for ¾" (19mm) WHITMEYER Style Ball Fitting	Stabiliser Mounting Bracket and Screws for Comfort Company Style Fitting



ADJUSTABLE FOREHEAD STABILISER



THE FOREHEAD STABILISER CAN BE FITTED TO POWERED WHEELCHAIRS, MANUAL WHEELCHAIRS AND SHOWER CHAIRS



**Fully
Adjustable in
Length, Angle
and Rotation**

A close-up, side-profile view of a blue exercise bike. The bike features a large, black, oval-shaped seat mounted on a black adjustable post. The main frame is a vibrant blue. A black pedal is visible in the foreground, and a black foot pedal with a strap is attached to the crank arm. A black foot pedal with a strap is also visible. A black foot pedal with a strap is also visible.

**Good
Positioning
and Comfort
Combined**

Leg and Foot Supports

Gel Ovation's leg, foot and knee protectors come in a variety of configurations and sizes to interface with all major wheelchair brands.

These products are designed to offer the ultimate in comfort and protection to knees, legs, heels and soles.

Our Gel Angle Adjustable Foot buckets

can be easily adjusted to the specific position and angle required. The padding is made from our unique Dimensional Gel for added comfort and protection against pressure issues.

Our Contoured Calf Support trough's padding

is made from our unique Dimensional Gel for added comfort and protection against pressure issues. This product comes with a removable cover and fixing bolts.



GEL FOOT PROTECTORS



GFPCA	GFPC	GFPCS
Gel & Foam Foot Protector Covers with Adhesive Back	Gel Foot Protectors Wooden Back with T-Nuts	Gel Single Piece Footplate Protector Wooden Back with T-Nuts
8" x 6" x 1.5" (200mm x 150mm x 38mm)	8" x 6" x 1.5" (200mm x 150mm x 38mm)	11" x 9" x 1.5" (280mm x 230mm x 38mm)

UNIT: Pair

UNIT: Pair

UNIT: Each

GEL KNEE PROTECTOR BUTTONS FOR ELEVATING LEG RESTS

KBE	KBET
Gel Knee Protector Buttons - Elevating Leg Rests	Gel Knee Protector Buttons - THIN Elevating Leg Rests
4" (100mm Diameter) 35mm Thick	4" (100mm Diameter) - 25mm Thick
Supplied with Pivot Bolts and Spacers	Supplied with Pivot Bolts and Spacers

UNIT: Pair

UNIT: Pair



GEL KNEE PROTECTOR BUTTONS



KBM	KBMT
Gel Knee Protector Buttons - Swing Away Leg Rests	Gel Knee Protector Buttons - THIN - Swing Away Leg Rests
4" (100mm Diameter) 35mm Thick	4" (100mm Diameter) 25mm Thick
Supplied with Tube Clamps and Bolts	Supplied with Tube Clamps and Bolts

UNIT: Pair

UNIT: Pair



KBMB
Gel Bariatric Knee Protector Buttons
8" x 6" x 1.5" (180mm x 150mm x 38mm) - Great for the Rea Assist Elevating Leg Rests

UNIT: Pair

GEL LEG PROTECTOR WRAPS



LLPW	LPW	LHC
Large Gel Leg Protector Wraps with Cut Out	Gel Leg Protector Wraps	Gel Leg Hanger Wrap
10.5" (270mm)	9.5" (240mm)	5.5" (140mm)

UNIT: Pair

UNIT: Pair

UNIT: Pair

GEL KNEE/SHOULDER PROTECTOR WRAPS

KPRT	KPR	KPRQ
Gel Knee/Shoulder Protector Wraps - Thin	Gel Knee/Shoulder Protector Wraps - Thick	Gel Knee/Shoulder Protector Wrap for larger than 1" Diameter Tubes & Corners
3.5" (90mm) Length 3/8" (8mm) Thick	3.5" (90mm) Length 5/8" (15mm) Thick	4.5" (115mm) Length 3/4" (20mm) Thick
UNIT: Pair	UNIT: Pair	UNIT: Pair



FULL LENGTH LEG PROTECTORS

FLGP10	FLGP12	FLGP14
Full-Length Leg Protectors 10" (260mm)	Full-Length Leg Protectors 12" (300mm)	Full-Length Leg Protectors 14" (360mm)
Supplied with Tube Clamps	Supplied with Tube Clamps	Supplied with Tube Clamps
UNIT: Pair	UNIT: Pair	UNIT: Pair



GEL ANGLE ADJUSTABLE FOOT BUCKET

AFBHO
Gel Angle Adjustable Foot Bucket
7" (180mm) High x 7" Deep and 6.5" (165mm) Wide
UNIT: EACH

DGCT
Contoured Calf Support Trough Padding
9 1/4" x 6 1/4" (235 x 155mm)
UNIT: PAIR



Our Gel Angle Adjustable Foot buckets can be easily adjusted to the specific position and angle required. The padding is made from our unique Dimensional Gel for added comfort and protection against pressure issues.

Our Contoured Calf Support trough's padding is made from our unique Dimensional Gel for added comfort and protection against pressure issues. This product comes with a removable cover and fixing bolts.

GEL CALF SUPPORT PANELS

GHS	GHSB
Gel Calf Support Panel with positioning Strap	Gel Calf Support Panel with positioning Strap
(Fits 16-20" (400 x 500mm) Wheelchair)	(Fits 20-24" (500 x 610mm) Wheelchair)

UNIT: Each

UNIT: Each



GHSBX	GHSP
Gel Calf Support Panel with positioning Strap	Gel Calf Support Panel with positioning Strap
(Fits 24-28" (610 x 710mm) Wheelchair)	(Fits 10-16" (260 x 400mm) Wheelchair)

UNIT: Each

UNIT: Each



GEL SPLIT CALF SUPPORT PANELS

GHS-2	GHSP-2
Gel Split Calf Support Panels with positioning straps	Gel Split Calf Support Panels with positioning straps
(Fits 16-20" (400 x 500mm) Wheelchair)	(Fits 10-16" (260 x 400mm) Wheelchair)

All UNITS: Each

All UNITS: Each



GHSP-2	GHSB-2
Gel Split Calf Support Panels with positioning straps	Gel Split Calf Support Panels with positioning straps
(Fits 20-24" (500 x 610mm) Wheelchair)	(Fits 24-28" (610 x 710mm) Wheelchair)

All UNITS: Each

All UNITS: Each



GEL FIXED KNEE ADDUCTORS WITH HARDWARE

These pads have the comfort of Gel and the strength of our hardware. These easy to fit pads can be placed in any position required and they STAY in position due to the additional set screws that come as standard.



SRABC	SRABD
Fixed Adductor Bracket with Gel Adductor	Fixed Adductor Bracket with Gel Adductor
Pad for 7/8" (22mm) Tubing	Pad for 1" (25mm) Tubing

ALL UNITS: Pair

ALL UNITS: Pair



SRABF
Fixed Adductor Bracket with Gel Adductor with Flat Mount Hardware for Seat Rail Systems. Compatible with Invacare, Sunrise Medical, Pride and Permobil

ALL UNITS: Pair

GEL SWING AWAY KNEE ADDUCTORS WITH HARDWARE



SRASC	SRASD
Swing Away or Removable Adductor Bracket with Gel Adductor	Swing Away or Removable Adductor Bracket with Gel Adductor
Pad for 7/8" (22mm) Tubing	Pad for 1" (25mm) Tubing

All UNITS: Pair

All UNITS: Pair



SRASF
As per our Swing away but with Flat Mount Hardware for Seat Rail Systems. Compatible with Invacare, Sunrise Medical, Pride and Permobil

All UNITS: Pair

STUMP SUPPORTS WITH HARDWARE

Our Stump Support Pads are slightly contoured and come with our unique Dimensional Gel. The ultimate protection against shear and heat build-up

67SRSSC	67SRSSD
Swing Away or Removable Stump Support Hardware for 7/8" (22mm) Tubing	Swing Away or Removable Stump Support Hardware for 1" (25mm) Tubing
Supplied with 6 x 7.5" (67DSSP) Dimensional Stump Support Pad	Supplied with 6 x 7.5" (67DSSP) Dimensional Stump Support Pad
All UNITS: Each	All UNITS: Each

89SRSSC	89SRSSD
Swing Away or Removable Stump Support Hardware for 7/8" (22mm) Tubing	Swing Away or Removable Stump Support Hardware for 1" (25mm) Tubing
Supplied with 8 x 9" (89DSSP) Dimensional Stump Support Pad	Supplied with 8 x 9" (89DSSP) Dimensional Stump Support Pad
All UNITS: Each	All UNITS: Each

67SRSSF	89SRSSF
As per swing away version but with Flat Mount Hardware for Seat Rail Systems. Compatible with Invacare, Sunrise Medical, Pride and Permobil	As per swing away version but with Flat Mount Hardware for Seat Rail Systems. Compatible with Invacare, Sunrise Medical, Pride and Permobil
Supplied with 6 x 7.5" (67DSSP) Dimensional Stump Support Pad	Supplied with 8 x 9" (89DSSP) Dimensional Stump Support Pad
All UNITS: Each	All UNITS: Each



REPLACEMENT STUMP SUPPORT PADS

89DSSP	67DSSP
Curved Stump Pad Support with Dimensional Gel	Curved Stump Pad Support with Dimensional Gel
8x9" (200 x 230mm)	6x7" (150 x 180mm)
ALL UNITS: Each	ALL UNITS: Each



GEL HEEL LOOP COVER

HLC3
Gel Heel Loop Cover with 3 Gel Sections
5.5" (140mm) Long
UNIT: PAIR



GEL LEG STRAPS

GLS26	GLS32	GLS36
Gel Leg Straps	Gel Leg Straps	Gel Leg Straps
55" (1397mm) Gel Leg Strap with Gel 3.5" (89mm) x 26" (660mm)	55" (1397mm) Gel Leg Strap with Gel 4" (100mm) x 32" (812mm)	55" (1397mm) Gel Leg Strap with Gel 4" (100mm) x 36" (914mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each



MULTI-AXIS HIP GUIDES

Gel Ovations Multi-Axis Hip Guides are perfect for positioning Hips and Thighs just where you need with the added comfort of Gel. They have the ability to swing in or out and can also be angled helping create the perfect position. Great when clients require leg support with a central footplate.

Support and positioning of the upper leg is now a lot easier and more comfortable.

These 6 axis adjustable hip guides with aluminium cones mount easily to the side of virtually any seating system and can be easily adjusted up and down and can be locked in the position of choice.

They are easy to remove and remain in the correct position.

Stainless steel on the primary support parts means they won't bend under tone. We've also added two (2) interlocking crown pivots to allow adjustment in and out. The track on the pad allows continuous forward and back adjustment.

**Help support
and position the
upper leg**





**Easy Six Axis
Positioning**

HG310	HG312	HG314	HG316
Multi-Axis Hip Guide	Multi-Axis Hip Guide	Multi-Axis Hip Guide	Multi-Axis Hip Guide
3.5 x 10" (90 x 260mm)	3.5 x 12" (90 x 300mm)	3.5 x 14" (90 x 360mm)	3.5 x 16" (90 x 400mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each

MOUNTING BRACKET FOR HIP GUIDE

Compatible with all seat rail systems including Sunrise Medical, Invacare, Pride and Permobil, when ordered with the HGMB Mounting Bracket. Please note that for the Salsa Seat System the Mounting Bracket is not required.

HGMB
Mounting Bracket for Hip Guide
ALL UNITS: Each



**Gel Ovation
Arm and Elbow
Supports**



Arm and Elbow Supports

These Gel Arm Pads offer the great comfort of Gel with a tough Neoprene cover for long-lasting life and durability. These pads are supplied with various M6 thread fixings that fit a variety of chairs.

Visit **www.gelovationseurope.com** for hole spacing

***“Our G Series Arm Pads
are soft, comfortable and
highly durable”***

G Series

GEL ARM REST PADS

310G	312G
Gel Armrest Pads 3.5 x 10" (90 x 260mm)	Gel Armrest Pads 3.5 x 12" (90 x 300mm)
c/w M6 T-Nuts	c/w M6 T-Nuts

ALL UNITS: Pair

ALL UNITS: Pair

314G	316G
Gel Armrest Pads 3.5 x 14" (90 x 360mm)	Gel Armrest Pads 3.5 x 16" (90 x 400mm)
c/w M6 T-Nuts	c/w M6 T-Nuts

ALL UNITS: Pair

ALL UNITS: Pair



210G	212G	214G
Gel Armrest Pads	Gel Armrest Pads	Gel Armrest Pads
2 x 10" (50 x 260mm) c/w M6 T-Nuts	2 x 12" (50 x 300mm) c/w M6 T-Nuts	2 x 14" (50 x 360mm) c/w M6 T-Nuts

ALL UNITS: Pair

ALL UNITS: Pair

ALL UNITS: Pair

STRENGTHENING PLATES FOR 314G & 316G

If you have clients who require a little extra strength in their arm pads these are the perfect solution. Our strengthening plates fit the 14" (360mm) and 16" (400mm) G Series Arm pads. They are particularly useful when you have stand transfer clients who exert extra pressure particularly on the ends of the pads.

314GBO	316GBO
Strengthening plate for 314G Arm Pads	Strengthening plate for 316G Arm Pads
ALL UNITS: Each	ALL UNITS: Each



TMS SERIES ARM PADS - FOR TOTAL ADJUSTMENT AND EASE OF FIT

These pads offer all the comfort of the G Series Pads but are also supplied with Labac Style tracking. The tracking is supplied with sliding T-Nuts and three lengths of bolt. This offers multiple fixing options with total ease of adjustment.



310TMS	312TMS	314TMS	316TMS
Gel Armrest Pads 3.5 x 10" (90 x 260mm)	Gel Armrest Pads 3.5 x 14" (90 x 300mm)	Gel Armrest Pads 3.5 x 16" (90 x 400mm)	Gel Armrest Pads 3.5 x 16" (90 x 400mm)
c/w- LaBac Style Mounting	- c/w- La Bac Style Mounting	- c/w- LaBac Style Mounting	- c/w- LaBac Style Mounting
ALL UNITS: Pair	ALL UNITS: Pair	ALL UNITS: Pair	ALL UNITS: Pair

210TMS	212TMS	214TMS
Gel Armrest Pads 2 x 10" (50 x 260mm)	Gel Armrest Pads 2 x 12" (90 x 300mm)	Gel Armrest Pads 2 x 14" (50 x 360mm)
- c/w- LaBac Style Mounting	- c/w- LaBac Style Mounting	- c/w- LaBac Style Mounting
ALL UNITS: Pair	ALL UNITS: Pair	ALL UNITS: Pair

8TMS	13TMS
LaBac Style Track with T-Nuts and Screws 8.5" (215mm)	LaBac Style Track with T-Nuts and Screws 13" (330mm)
ALL UNITS: Each	ALL UNITS: Each



GT SERIES - OUR TOTALLY MODULAR ARM PAD SYSTEM

These Gel Arm Pads are modular offering your client a totally personal service. This allows you to create the support and comfort that they require. The accessories include waterfall pads, trough pads, and elbow pads. The GT Series arm pads also come with five rows of Labac Style Tracking giving you more adjustment than ever before.

Please choose from four different arm pad sizes and combine with the required accessories.

1. First choose the arm pad size that you require:



**GT Series
Arm Pad with
Labac Style Tracking
that allows ease of
adjustment and option
to fit accessories.**

408GT	410GT	412GT	415GT
Flat Gel Pad	Flat Gel Pad	Flat Gel Pad	Flat Gel Pad
c/w Removable Cover & LaBac style Tracking 4 x 8" 100 x 200mm	c/w Removable Cover & LaBac style Tracking 4 x 10" 100 x 260mm	c/w Removable Cover & LaBac style Tracking 4 x 12" 100 x 300mm	c/w Removable Cover & LaBac style Tracking 4 x 15" 100 x 380mm
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each

2. Next choose the required accessories:

PF182	PF185	PF187
Modular Trough Pad 6" (150mm)	Modular Trough Pad 8" (200mm)	Modular Trough Pad 10" (260mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each



PF183	PF186	PF188
Modular Waterfall Pad 6" (150mm)	Modular Waterfall Pad 8" (200mm)	Modular Waterfall Pad 10" (260mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each



PF210	PF211
Contoured Elbow Stop Kit	Contoured Elbow Stop Kit
Left Hand	Right Hand
ALL UNITS: Each	ALL UNITS: Each



Please note that these Elbow Stop Kits wrap around the corner of the Arm Pads by 1.5" (40mm). Important to remember when ordering a Trough Pad in combination with an Elbow Stop.

PF215	PFMPS	PFMPL
GT Series Hand Pad	Small Mouse Pad	Large Mouse Pad
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each





GWT SERIES ARM PADS - THE WATERFALL RANGE

These heavy-duty arm pads are perfect for power chairs where clients require that little extra padding. These pads are also supplied with the Labac Style Tracking so you can be assured it will fit almost any chair. The drop of the arm is 3" (75mm).

410GWT	412GWT
Waterfall Style Gel Pad with LaBac style Mounting	Waterfall Style Gel Pad with LaBac style Mounting
4 x 10" (100 x 260mm)	4 x 12" (100 x 300mm)
ALL UNITS: Each	ALL UNITS: Each

414GWT	416GWT
Waterfall Style Gel Pad with LaBac style Mounting	Waterfall Style Gel Pad with LaBac style Mounting
4 x 14" (100 x 360mm)	4 x 16" (100 x 400mm)
ALL UNITS: Each	ALL UNITS: Each



UNIVERSAL ELBOW STOP BRACKET WITH FLAT GEL PAD

These Universal Gel Elbow Stop Brackets can be fitted to our G -Series (310G, 312G ,314G or 316G), GT Series or GWT Series Arm Pads.

35UES	35UES
45 Degree elbow stop bracket with curved Dimensional Gel Pad	45 Degree elbow stop bracket with curved Dimensional Gel Pad
3" x 5" (100mm x 130mm)	4" x 6" (100mm x 150mm)

ALL UNITS: Each

ALL UNITS: Each

57UES	UESF
45 Degree elbow stop bracket with curved Dimensional Gel Pad	45 Degree Elbow Stop Bracket with flat Gel Pad
5" x 7" (130mm x 180mmn)	3" x 4" (75mm x 100mm)

ALL UNITS: Each

ALL UNITS: Each

UESC
45 Degree Elbow Stop Bracket with CURVED Gel Pad
3" x 4' (75mm x 100mm)

ALL UNITS: Each



BARIATRIC UNIVERSAL ELBOW STOP BRACKETS

These heavy-duty elbow stop brackets come with a 30-degree bend that will not bend out of shape in normal use, particularly useful for clients with high tone. Compatible with all Gel Ovation elbow stop pads from 3 x 4" (75 x 100mm) to 5 x 7" (130 x 180mm).

HDESUL	HDESUS
Heavy Duty Elbow Stop Hardware	Heavy Duty Elbow Stop Hardware
for use with 10" - 12" (260mm - 300mm) Arm Pads	for use with 14" - 16" (360mm - 400mm) Arm Pads
ALL UNITS: Each	ALL UNITS: Each



REPLACEMENT / GEL ELBOW STOPS

Order the BARIATRIC UNIVERSAL ELBOW STOP BRACKET and then choose what size Gel Elbow Stop you require.

35DGE1	46DGE1	57DGE1
Dimensional Gel Elbow Pad	Dimensional Gel Elbow Pad	Dimensional Gel Elbow Pad
3 x 5" (75 x 130mm)	4 x 6" (100 x 150mm)	5 x 7" (130 x 180mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each

67DGE1	89DGE1
Dimensional Gel Elbow Pad	Dimensional Gel Elbow Pad
6 x 7.5" (150 x 190mm)	8 x 9" (200 x 230mm)
ALL UNITS: Each	ALL UNITS: Each



GEL OVATIONS FIXED POSITION ELBOW STOP BRACKETS

EBLSD	EBLSC
Fixed position elbow stop with 1" (25mm) round tube adaptor Gel pad	Fixed position elbow stop with 7/8" (22mm) round tube adaptor Gel pad
3" x 4" (80mm x 100mm)	3" x 4" (80mm x 100mm)

ALL UNITS: Pair

ALL UNITS: Pair



GEL OVATIONS ADJUSTABLE ELBOW STOP BRACKETS

SRPESD	SRPESC
Position adjustable elbow stop for 1" (25mm) tubing Gel pad	Position adjustable elbow stop for 7/8" (22mm) tubing Gel pad
3 x 4" (80 x 100mm)	3 x 4" (80 x 100mm)

ALL UNITS: Pair

ALL UNITS: Pair



SRES	SRESC
Position adjustable swing away or removable elbow stop for 1" (25mm) tubing. Gel pad	Position adjustable swing away or removable elbow stop for 7/8" (22mm) tubing Gel pad
3" x 4" (80 x 100mm)	3x 4" (80 x 100mm)

ALL UNITS: Pair

ALL UNITS: Pair



EXTENSION STEMS FOR ELBOW STOP BRACKETS

EES2	EES3	EES24
2" (50mm) extension stem	3" (75mm) extension stem	4" (100mm) extension stem
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each



'Gel Ovation's premium hardware includes serrations, set screws and interlocking crowns or cones. This ensures that the hardware will not slip, bend or break which protects the users comfort and postural needs'



Lateral Supports
and Hardware



Lateral Supports and Hardware

There are three steps to follow to order your Lateral Pad Solution:

- 1) Choose the Size of Lateral Pad required.
- 2) Decide on a Left or Right Lateral Support Hardware
- 3) Order the correct size tube clamp (P44) to attach the lateral support hardware to the back support tubes of the wheelchair

- Manufactured from Medical Grade Silicone
- Unique Design Reduces Shear & Friction Forces
- Keeps Cool
- Lateral Pads available from 3 x 5" (75 x 130mm) up to 8 x 9" (200 x 230mm)



DIMENSIONAL GEL LATERAL PADS

35DGLP1	46DGLP1
Dimensional Gel Lateral Pads	Dimensional Gel Lateral Pads
3 x 5" (75 x 130mm)	4 x 6" (100 x 150mm)

UNIT: Each

UNIT: Each



57DGLP1	67DCLP1	89DCLP1
Dimensional Gel Lateral Pads	Dimensional Gel Lateral Pads	Dimensional Gel Lateral Pads
5 x 7" (130 x 180mm)	6 x 7.5" (150 x 190mm)	8 x 9" (200 x 230mm)

UNIT: Each

UNIT: Each

UNIT: Each



ARTICULATING LATERAL SUPPORT HARDWARE

SRRLHL	SRRLHR
Articulating Lateral support hardware	Articulating Lateral support hardware
Left Hand	Right Hand

UNIT: Each

UNIT: Each



Simply Lift Bracket to Swing Away



TUBE CLAMPS TO FIT LATERAL HARDWARE AND OTHER ACCESSORIES



SAB2D1	SAB2C1	SAB3D1	SAB3C1
Gel Ovations accessory clamp	Gel Ovations accessory clamp	Gel Ovations accessory clamp	Gel Ovations accessory clamp
2" long for 1" (25mm) round tubing	2" long for 7/8" (22mm) round tubing	3" long for 1" (25mm) round tubing	3" long for 7/8" (22mm) round tubing

ALL UNITS: Each

ALL UNITS: Each

ALL UNITS: Each

ALL UNITS: Each



TC1	TC2
1" (25mm) Serrated Edge Tube Clamp	2" (50mm) Serrated Edge Tube Clamp
1" (25mm)	2" (50mm)

ALL UNITS: Each

ALL UNITS: Each



SR75ADP
These sleeves allow you to fit 7/8" (22mm) Clamps on to 3/4" (19mm) Tubing
Tube Extension Sleeve for 3/4" (19mm) Tubing

ALL UNITS: Set

**Wraps, Covers
and Stick-on
Pads**



Wraps, Covers and Stick-on Pads

Metal framing, seat belts and straps all have one thing in common; they offer a hard, unyielding surface against which skin and soft tissue have little or no chance of surviving intact.

Gel Ovation make a variety of products designed to offer the client comfort and protection whenever and wherever contact is inevitable by providing a variety of attachment methods. We are able to offer a pressure-relieving surface that can eliminate problem surfaces - from seat belts to metal tubing.

- Easy to fit
- Instant Solutions
- Blends in with wheelchair or mobility aid
- Simply Velcro or stick into position
- Bespoke service available



GEL PROTECTOR TUBE WRAPS

PTW6	PTW8
Gel Protector Tube Wraps	Gel Protector Tube Wraps
6" (150mm)	8" (200mm)
ALL UNITS: Pair	ALL UNITS: Pair



GEL DRIVE HANDLE COVER

THC
Gel Drive Handle Cover
4" (100mm)
UNITS: Each



GEL STRAP COVER

GSC5	GSC8
Gel Strap Cover with Velcro Closure + Grommets	Gel Strap Cover with Velcro Closure + Grommets
5" (130mm)	8" (200mm)
ALL UNITS: Pair	ALL UNITS: Pair



GEL SEAT BELT COVERS

SBCV	SBCPV	SBCBV
Gel Seat Belt Cover Standard Adult	Gel Seat Belt Cover Paediatric	Gel Seat Belt Cover Bariatric
ALL UNITS: Set	ALL UNITS: Set	ALL UNITS: Set



BACK CANE COVER



BCC16	BCC18	BCC20
Back Cane Cover	Back Cane Cover	Back Cane Cover
16" (400mm)	18" (460mm)	20" (500mm)
UNIT: Each	UNIT: Each	UNIT: Each

GEL PADDED SLIP OVER BACK REST COVER

These Gel Slip Over back covers are compatible with almost all makes and model of wheelchair. They simply slide over the push handles and the back is secured into position by using the Velcro straps at the rear.

GPSB16	GPSB18	GPSB20
Gel Padded Slip over Back Cover 16" (400mm)	Gel Padded Slip over Back Cover 18" (460mm)	Gel Padded Slip over Back Cover 20" (500mm)
ALL UNITS: Each	ALL UNITS: Each	ALL UNITS: Each



GEL WALKER, ROLLATOR AND CRUTCH HANDLE COVERS

WHC4	WHC5
Gel Walker / Crutch Handle Cover - 4" (100mm)	Gel Walker / Crutch Handle Cover - 5" (130mm)
ALL UNITS: Pair	UNIT: Pair



RHC
Gel Rollator Handle Cover
UNIT: Pair



STICK-ON PADS

Arms, legs and knees don't just come in contact with parts of the wheelchair. Desktops, vehicle surfaces and standing frames all provide an opportunity for skin abrasions and pressure sores. Stick-on pads, while offering all of the comfort and safety of all our gel products, also provide the convenience of being easily applied to any flat surface, whenever and wherever needed.



ESP	SESP	34SP	48SP
Egg-shaped Stick-on Pad	Small Egg-shaped Stick-on Pad	Gel Stick-on Pads 3 x 4" (80 x 100mm)	Gel Stick-on Pads 4 x 8" (100 x 200mm)
UNITS: Each	UNITS: Each	UNIT: Pair	UNITS: Pair

412SP	1575SP	2075SP	2575SP
Gel Stick-on Pads 4.5" x 12" (120 x 300mm)	Flat Stick-on Pads with Rounded End 1.5" x 7.5" (40mm x 190mm)	Flat Stick-on Pads with Rounded End 2" x 7.5" (50mm x 190mm)	Flat Stick-on Pads with Rounded End 2.5" x 7.5" (65mm x 190mm)
UNIT: Pair	ALL UNITS: Pair	ALL UNITS: Pair	ALL UNITS: Pair

**Gel Ovation
Shower Chair
and Bathroom
Accessories**



Gel Ovation Aqua Range of Posture and Pressure Care Solutions

Gel Ovation provide many simple Gel Wraps and Covers for use in the wet room environment. This helps to protect skin as part of a 24-hour care programme when it is at its most vulnerable.

SHOWER CHAIR ACCESSORIES



LPWA	KPRA	TACA	INVSCFP
Gel Leg Protector Wraps	Gel Shoulder / Knee Protector Wraps	Gel Tubular Arm Protectors	Gel Ocean Shower Chair Foot Protectors with Velcro Straps

ALL UNITS: Pair

ALL UNITS: Pair

UNIT: Pair

UNIT: Pair



INVSCFPIP	INVSCLP	INVCHR	INVSCBR
Gel Ocean Shower Chair One Piece Foot Protector with velcro straps	Gel Ocean Shower Chair Gel Lateral Protectors with Velcro Straps	Gel Ocean Shower Chair Gel Headrest Wrap	Gel Ocean Shower Chair Soft Neoprene Back Rest

UNIT: Each

UNIT: Pair

UNIT: Each

UNIT: Each

GEL OVATION TOILET SEAT AND SHOWER CHAIR OVERLAYS

Gel Ovation manufacture many bespoke toilet seat and shower chair seat overlays. Using our Dimensional Gel, with a VP Cover and Non Slip Base, we offer a low profile Pressure Solution for your personal care needs. For more information regarding our bespoke service please visit: www.gelovationseurope.com



OSCSO

Ocean Shower Chair Seat Overlay

UNIT: Each



GTSO

Geberit Toilet Seat Overlay

UNIT: Each

GTSOM


Geberit Meyra Toilet Seat Overlay

UNIT: Each

AQUA HARNESSES EXTRA SOFT FOR COMFORT

Designed for people with trunk instability, Chest Support Harnesses are made of elastic neoprene that is kind to the skin and provides greater freedom of movement.

Each adjustment of the fixing straps results in perfect adaptation to the clinical requirements. The Chest Support Harnesses are easy-care and machine-washable up to 30° C and are available in three sizes in the standard version.



SIZING CHART FOR THE AQUA HARNESSES AND HIP BELTS

PRODUCT	CODE	A	B	C	D
Aqua Chest Support Harness SMALL	AHS	33.0 cm	22.0 cm	13.0 cm	20.5 cm
Aqua Chest Support Harness MEDIUM	AHM	38.0 cm	25.0 cm	15.0 cm	23.0 cm
Aqua Chest Support Harness LARGE	AHL	42.0 cm	27.0 cm	17.0 cm	24.5 cm

LAP BELTS

Gel Ovation provides Soft Neoprene, 2 Point and 4 Point, Lap Belts for positioning requirements whilst showering. These padded belts help to protect the skin and offer support when it is at its most vulnerable.

Please note that these belts are not safety belts and should not be used for transportation.

2-POINT NEOPRENE LAP BELT WITH PLASTIC BUCKLE



A2PBS	A2PBM	A2PBL
Aqua 2-point Neoprene Lap Belt with Plastic Buckle	Aqua 2-point Neoprene Lap Belt with Plastic Buckle	Aqua 2-point Neoprene Lap Belt with Plastic Buckle

SMALL

MEDIUM

LARGE

Please visit www.gelovationseurope.com for full dimensions

4-POINT NEOPRENE LAP BELT WITH PLASTIC BUCKLE



A4PBS	A4PBM	A4PBL
Aqua 4-Point Neoprene Lap Belt with Plastic Buckle	Aqua 4-Point Neoprene Lap Belt with Plastic Buckle	Aqua 4-Point Neoprene Lap Belt with Plastic Buckle

SMALL

MEDIUM

LARGE

Please visit www.gelovationseurope.com for full dimensions

GEL OVATIONS FOREHEAD STABILISER FOR THE OCEAN SHOWER CHAIR

Our Forehead Stabiliser can be fitted to the Ocean Shower Chair by adding a simple Adaptor. This gives multiple fixing positions to accommodate the users needs during personal care. For further information on the Forehead Stabiliser please see Page 14.



SCHA	SRAFS510	SRAFS512	SRAFS514
Ocean Shower chair Headrest Adapter - Ball Joint for extra adjustment during fitting (Must be used in conjunction with the Ocean Ergonomic Head Rest Holder)	Adjustable Forehead Stabiliser with 5" x 10" (130 x 260mm) Headrest Pad	Adjustable Forehead Stabiliser with 5" x 12" (130 x 305mm) Headrest Pad	Adjustable Forehead Stabiliser with 5" x 14" (130 x 360mm) Headrest Pad
UNIT: Each	UNIT: Each	UNIT: Each	UNIT: Each







Gel Ovation Europe Ltd
Charlton Road
Bristol BS10 6NF
United Kingdom

Tel: +44 (0)1454 285 071
Email: info@gelovationseurope.com
Visit: www.gelovationseurope.com



Copyright 2019 All contents copyright Gel Ovation Europe Ltd. All copyright, design rights and other intellectual property in Gel Ovation products and designs are and remain the property of Gel Ovation Europe Ltd. All product details correct at time of going to print. Due to continual product development we reserve the right to change specifications or materials without prior notice.

VAT Number: 921 3722 49
Company Number: 6362508

Issue Dec 2019