



# GLS5.1

## Active lifter

*– an easy, comfortable and safe way to raise and support patients with partial or complete control of upper body*

### **One lift for many functions**

- Sit to stand
- Standing balance (with foot tray removed)
- Transfers from chair to chair
- Transfers from bed to chair
- Transfers from bed to toilet

### **Adjusts easily to comfortably fit most heights and sizes**

- Padded and adjustable knee support
- Multiple sling attachment points
- Ergonomically profiled push handle for all caregiver sizes
- Multiple sling sizes to create optimum comfort and support (option)
- Add-on hip sling to assist with hips and buttocks raising (option)

### **Easy and safe to use and understand**

- Electric up, down, and leg width adjustment
- Dual operations options – hand control or control box
- Rear foot push bar to increase maneuverability
- Manually and electrically operated emergency lowering
- Easily accessible emergency stop



### Strong, reliable and powerful

- SWL 155 or 205 kg
- Designed for intensive use in professional environments
- Meets and exceeds all current norms and regulations in industry
- Strong 25,9 V/2,25 AH Li-Ion battery ensures long usage time, even in heavy use environments
  - Shorter charging time to full capacity
  - Prolonged battery life
  - More lifts from a single charge
- Strong, noiseless and fast Linak actuators
- Varied charging options – on board lift and off board lift via extra battery
  - Built in and easily accessible docking station for battery
  - Extra battery packs as option
  - Fully automatic switch mode charger
- Removable and easy to clean foot tray
- The battery capacity is indicated in 3 intervals in the hand control;
  - Low 0 - 25%
  - Medium 25 - 50%
  - High 50 - 100%Audio signal at low battery voltage when approx. 25%.
- Visual service indicator lights up when maintenance is required.

Time to care

V. Guldmann A/S  
Tel. +45 8741 3100  
info@guldmann.com  
www.guldmann.com