







A photograph showing a person with white hair, wearing a yellow and green striped shirt, seated in a wheelchair. A hand is visible on the left, holding a blue and black SANO People Carrier device, which is being used to assist the person. The device is attached to the wheelchair's frame. The background is a plain, light-colored wall.

**Quick and easy to the correct
People carriers**

	LIFTKAR [®] PTR		LIFTKAR [®] PT	
	PTR	PTR-L	PT S	PT Outdoor
				
Description	easy-operate powered tracked stair climber	longer version (+132 mm) of the easy-operate powered tracked stair climber	safety and ergonomic handling for indoors	safe transport outside the home, too
Ideal for	universal wheelchairs	lightweight electric wheelchairs	climbs straight staircases (up and down) indoors	uneven terrain outdoors
Load capacity	130 kg / 160 kg	130 kg / 160 kg	130 kg / 160 kg	120 kg / 150 kg
Total weight	38.2 kg	39.4 kg	35 kg	37 kg
Weight of the heaviest part	28.6 kg (base unit)	29.8 kg (base unit)	21.5 kg (climber unit with seat)	24 kg (climber unit with seat)
Step height / angle (max.)	35°	35°	220 mm	220 mm
Lightweight construction	✓	✓	✓	✓
Operates on steep and narrow spiral staircases	–	–	✓	✓
Easy to dismantle and stow	✓	✓	✓	✓
Integrated seat	–	–	✓	✓
Wheelchair transporters	✓	✓	–	–
Single-step mode & continuous mode	–	–	✓	✓
Acoustic and visual inclination warning	✓	✓	✓	✓
Motor-supported assistance for getting up	–	–	✓	✓
Reliable step edge brakes	–	–	✓	✓
Can be used on a wide range of wheelchairs	✓	✓	–	–
Autonomy (max.)	up to 1,000 Stufen	up to 1,000 Stufen	up to 500 steps	up to 500 steps
Climbing speed	4.5 / 5.5 / 6.5 (m/min)	4.5 / 5.5 / 6.5 (m/min)	10 / 14 / 18 (steps/min)	10 / 14 / 18 (steps/min)
Landing dimensions (min.)	970 mm x 970 mm	1,150 mm x 1,150 mm	800 mm x 900 mm	800 mm x 900 mm

	PT Fold	PT Uni	PT Adapt	PT Plus
				
Description	the ideal solution for patient transport	the universal solution for hand-propelled and transport wheelchairs	to fit to the existing wheelchair	the complete solution with a high-quality wheelchair
Ideal for	patient transport and very narrow spiral staircases	professional transport of people with a variety of wheelchairs	adapter ensures perfect clamping of the existing wheelchair	simple and easy docking of the wheelchair provided
Load capacity	130 kg / 160 kg	130 kg / 160 kg	130 kg / 160 kg	125 kg
Total weight	28 kg	28 kg	23 kg	22.5 kg (without wheelchair)
Weight of the heaviest part	–	18 kg (climber with platform)	14.5 kg (climber)	14.5 kg (climber)
Step height / angle (max.)	220 mm	220 mm	220 mm	220 mm
Lightweight construction	✓	✓	✓	✓
Operates on steep and narrow spiral staircases	✓	✓	✓	✓
Easy to dismantle and stow	✓	✓	✓	✓
Integrated seat	✓	–	–	–
Wheelchair transporters	–	✓	✓	✓
Single-step mode & continuous mode	✓	✓	✓	✓
Acoustic and visual inclination warning	✓	✓	✓	✓
Motor-supported assistance for getting up	✓	✓	✓	✓
Reliable step edge brakes	✓	✓	✓	✓
Can be used on a wide range of wheelchairs	–	✓	✓	–
Autonomy (max.)	up to 500 steps	up to 500 steps	up to 500 steps	up to 500 steps
Climbing speed	10 / 14 / 18 (steps/min)	10 / 14 / 18 (steps/min)	10 / 14 / 18 (steps/min)	10 / 14 / 18 (steps/min)
Landing dimensions (min.)	800 mm x 900 mm	1,100 mm x 900 mm	900 mm x 900 mm	900 mm x 900 mm