



An effective aid in the prevention and treatment
of pressure ulcers ^(2, 4)

1 Repose® Mattress Overlay

Product order code: 6101100

Dimensions: L 1780mm x W 770mm

WL: 139kg



2 Repose® Mattress Overlay & Cushion set

Product order code: 6121100

Dimensions: L 1780mm x W 770mm
and L 450mm x W 450mm

WL: 139kg



3 Repose® Trolley Mattress Overlay

Product order code: 6601100

Dimensions: L 1780mm x W 550mm

WL: 139kg



4 Repose® Cushion

Product order code: 6201100

Dimensions: L 450mm x W 450mm

WL: 139kg



5 Repose® Foot Protector

Product order code: 6501100

Dimensions: L 380mm x W 200mm x H 180mm

WL: 139kg



6 Repose® Paediatric Mattress Overlay

Product order code: 6801100

Dimensions: L 1290mm x W 670mm

WL: 139kg



7 Repose® Babytherm

Product order code: 6401100

Dimensions: L 760mm x W 530mm

WL: 139kg



8 Repose® Babynest

Product order code: 6801101

Dimensions: L 580 mm x W 330 mm x H 100 mm

WL: 3kg



9 Repose® Wedge

Product order code: 6711100

Dimensions: L 690mm X W 450mm

WL: 139kg



10 Repose® Companion Cover

Product order code: 6990002

Dimensions: L 2050mm x W 880mm

WL: 139kg



11 Repose® Companion Sets

Single Set (x1)

Companion x 1
Mattress Overlay x 1
Pump x 1

Product order code:
6990003

Double Set (x2)

Companion x 2
Mattress Overlay x 2
Pump x 2

Product order code:
6990004



12 Repose® Thermal Companion

A universal mattress providing patient warming, effective pressure redistribution, and easy lateral transfer



Repose Thermal Carbon Pad

Dimensions: L 1800mm x W 750mm x H 20mm

Repose Thermal Companion Control Unit

Dimensions: L 285mm x W 125mm x H 150mm

Repose Companion Cover

Dimensions: L 2050mm x W 880mm

Repose Mattress & Pump

Dimensions: L1900mm x W 770mm

Repose® Total Care Packs

Acute

Product order code: 6882995

Rental code: 6882998

Mattress Overlay x 10
Cushion x 7
Foot Protector (Pair) x 5
Wedge x 1
Pump x 2

Nursing Home

Product order code: 6882997

Mattress Overlay x 5
Cushion x 3
Foot Protector (Pair) x 2
Pump x 1



Community Care Pack

Product order code: 6882996

Mattress Overlay x 5
Cushion x 3
Foot Protector (Pair) x 2
Pump x 5



Repose® Refresh Covers

Refresh Mattress Cover

Product order code: 6050101

Dimensions: L 1780mm x W 770mm

WL: 139kg

Refresh Cushion Cover

Product order code: 6231101

Dimensions: L450mm x W 450mm

WL: 139kg



Repose® Multi Stretch Moisture Vapour Permeable Cover

Trolley Mattress Cover

Product order code: 6980001

Dimensions: L 1880mm x W 580mm

WL: 139kg

Cushion Cover

Product order code: 6211100

Dimensions: L 520mm x W 520mm

WL: 139kg

An effective aid in the prevention and treatment of pressure ulcers

Indications for use

Repose products provide effective pressure redistribution for all persons at risk of developing pressure ulcers, including those assessed as very high risk. Repose is also appropriate for patients with pressure related tissue damage in-line with NPUAP-EPUAP Pressure Ulcer Treatment Clinical Practice Guidelines 2009 - clinical supervision is advised where the damage is severe. ^(1,2,8,13)

Limitations for use

- Repose is not suitable for the following:
- Persons weighing in excess of 139kg
 - Persons with unstable spinal fractures
 - Where body shape is such that the person cannot be fully supported by the Repose mattress overlay or cushion
 - The Repose mattress overlay is suitable for use on top of a bed or trolley mattress. Use side rails in accordance with manufacturer's instructions and appropriate risk assessment
 - The Repose cushion is intended for use on a static chair fitted with arm rests

Precautions

- If deterioration in skin condition is noted, clinical advice should be sought and, if so advised, discontinue use of product
- Avoid direct contact with sharp objects
- When deflated, keep away from babies and small children
- Keep away from pets

* All Repose products are latex free

REPOSE CAN BE CLEANED AND REUSED



REPOSE CAN BE CLEANED WITH SOAP OR DETERGENT AND WATER. REPOSE CAN BE DISINFECTED IN ACCORDANCE WITH INFECTION CONTROL POLICY FOR STANDARD MATTRESSES. ALCOHOL & PHENOL BASED DISINFECTANTS SHOULD NOT BE USED. ALLOW TO DRY BEFORE USE OR STORAGE

REPOSE IS SUPPLIED WITH A MANUFACTURERS WARRANTY FOR 12 MONTH FROM DATE OF DISPATCH

Repose® is clinically effective, cost effective, comfortable, durable and requires no maintenance.

Repose® products provide effective pressure redistribution for all persons at risk of developing pressure ulcers, including those assessed as very high risk. Repose is also appropriate for patients with pressure related tissue damage in-line with NPUAP-EPUAP Pressure Ulcer Treatment Clinical Practice Guidelines 2009 - clinical supervision is advised where the damage is severe. ^(1,4,5,6)

Clinically effective – comparable efficacy to significantly more expensive powered mattresses.

Authors reported in one RCT involving 80 elderly patients with fractured neck of femur no statistically significant difference was found between the Repose and Nimbus* ⁽¹⁾

Cost effective – up to 50 more patients benefit from Repose than conventional powered mattresses.

"There has been a definite shift away from the use of dynamic systems. Due to the increased availability of Repose mattresses, patients are being 'upgraded' more quickly. ... Although cost impact has not been looked at, the Trust was spending in excess of £300,000 per year on dynamic systems. The budget for all pressure relieving equipment is now £62,000 per annum" ⁽²⁾

Durable – average life span equivalent to more expensive dynamic systems

Two large community equipment loan stores report Repose mattress average life cycle to be in excess of 20 months. The new Refresh cover provides the opportunity to further extend the life of a Repose and maintain pristine condition ⁽³⁾

Portable – lightweight, dual purpose pump packaging – easy to store and transport.

Popular with NHS Supply Chain logistics, Community Equipment Services loan stores and Nurses

Easy to use – clear instructions, quick 'fool proof' inflation, easy to set up ready for use.

Patented 'smart valve' technology inside the pump ensures the product is primed and maintained at the optimum effective pressure.

Comfortable – a silent, unobtrusive system provides major positive benefits to the patient's quality of life.

Research by Price et al (2003) demonstrated that Repose aids patient comfort. *"In this pilot study statistically significant improvements in sleep and pain were noted over a four-week period"* ⁽⁷⁾



1. Price. P et al. 1999. Challenging the pressure sore paradigm. Journal of Wound Care, April, Vol 8, No 4.
2. Ballard Wilson, A. Pressure ulcer prevention audit: What are the benefits of doing it? Poster EPUAP Berlin, 2006.
3. Staffordshire County Council Loan Store 2007. Data on file.
4. Bale. S. Interim results of a randomised controlled trial of two pressure - relieving systems in the prevention of pressure sores in patients with fractured neck of femur. Presented at the 7th European Conference on Advances in Wound Management, Harrogate, 1997.
5. Osterbrink, J et al. Clinical evaluation of the effectiveness of a multimodal static pressure relieving device. 8th European Pressure Ulcer Advisory Panel Open Meeting, Aberdeen. May 5th – 7th, 2005.
6. NPUAP-EPUAP 2009 Pressure Ulcer Treatment QRG p19
7. Price P et al. 2003 The use of a new overlay mattress in patients with chronic pain: impact on sleep and self-reported pain. Clinical Rehabilitation 2003:17.

*Nimbus is a registered trademark of Huntleigh Technology PLC

Repose and Repose Companion are registered trade marks of Frontier Medical Group.
© September 2010 Frontier Therapeutics Limited. All rights reserved. GP0041_A20_05



Frontier Therapeutics Limited
Newbridge Road Industrial Estate
Blackwood, South Wales NP12 2YN

Tel: +44(0)1495 235800
Fax: +44(0)1495 235808
Email: info@frontiermedical.eu

www.frontiermedical.eu