Fellow Supreme











Suitable for many users, the Fellow Supreme combines high quality ergonomic design with a low centre of gravity for increased safety and stability. The forward facing handles allow the user to walk in a more upright position, which is better for both the back as well as the shoulders.

Unique Design

The Fellow Supreme features a unique fork construction that makes it extremely easy to maneuver, especially in smaller spaces. The flip-up seat allows the user to walk into the seat which allows them closer access to counters while doing work in and around the home.

Ergonomic Grips

The Fellow Supreme's soft braking system combined with the ergonomic design of the handles provide users with reduced hand functions the ability to easily lock the wheels in place and comfortably maneuver the rollator while walking. The Fellow Supreme folds easily with a one-hand operation and

can stand independently while folded. Its vast handle height adjustment range make it suitable for most users.

Practical Accessories

The Fellow Supreme comes with an assortment of accessories to make day-to-day life easier. From a tray to a cane holder, you can easily find an accessory to suit your needs. The grid easily removes for cleaning and maintenance.

HUMAN CARE SWEDEN (HQ)

Årstaängsvägen 21B 117 43 Stockholm Phone:- +46 8 665 35 00 Fax:- +46 8 665 35 10 info.se@humancaregroup.com

HUMAN CARE AUSTRALIA

Unit 4, 49 Bertie Street Port Melbourne, VIC 3207 Phone: +61 1300 132 736 Fax: +61 [03] 9645 4954 info.au@humancaregroup.com

HUMAN CARE CANADA

10-155 Colonnade Road Ottawa, ON K2E 7K1 Phone: +613.723.6734 Fax: +613.723.1058 info.ca@humancaregroup.com

HUMAN CARE NEDERLAND

Elspeterweg 124 8076 PA Vierhouten Phone: +31 577 412 171 Fax: +31 577 412 170 info.nl@humancaregroup.com

HUMAN CARE UNITED STATES

8006 Cameron Road, Suite K Austin, TX 78754 Phone: +512.476.7199 Fax: +512.476.7190 info.us@humancaregroup.com

Accessories/Features



TrayPractical for resting lightweight items.
Maximum load 2kg/4.4lb.
Part no. 10667



Mesh basket Ideal for those great shopping adventures. Maximum capacity 4.5kg/10lb. Part no.: 10666



Cane holder
Simplifies transport of a cane or umbrella.
Part no. 10603



Slow-Down Brake
Prevents the rollator from
moving too fast. Part
no. 10625 (set), 88121
[right], 88122 [left]



Soft Bag Soft carry case ideal for personal items. Part no. 08808



Light kitIlluminate your way while walking. Part no. 91091



Oxygen Tank Holder Part no. 10671



Backrest
Backrest for extra
comfort while seated.
Part no. 10668



Reflector Kit Includes 30 selfadhesive reflectors. Part no. 91090



Security Lock
Flex cord that protects
rollator from theft while
unattended.
Part no. 91092



One-Hand Brake Allows control of brakes with one hand. Must be ordered with rollator. Part no. 10669 (right), 10676 (left)



Newspaper Strap Keeps magazines and newspapers secure while walking. Part no. 10679

Technical Specifications

Part no.	85400
Handle Height (Min - Max)	73cm - 100cm / 28.7" - 39.4"
Seat Height	59cm/23.2"
Overall Width	63cm/24.8"
Overall Length	63.5cm/25"
Width between Handles	44cm/17.3"
Distance between Rear Wheels	55cm / 21.7"
Wheel Size	20cm x 4cm / 7.7" x 1.7"
Product Weight	8kg / 17.6lb
Maximum User Weight	150kg/330.7lb