# SOUND ASSIST CAN HELP IN COUNTLESS SITUATIONS

- Partner microphone mode
  One-to-one conversations with friends and family sound clear wherever you are.
- Table microphone mode Hear important meetings, lively dinners, and group conversations more easily.
- Hands-free phone calls Great sounding calls on the go.
- Streaming from any Bluetooth device Send music, film, or TV audio to your hearing aids.
- Remote control for Widex MOMENT Bluetooth hearing aids Control your sound, easily.
- Telecoil mode Telecoil functionality even if your hearing aids don't have telecoil.



SOUND LIKE NO OTHER



Widex is a registered trademark of WIDEX A/S. Printed by Rosendahls a/s. 2022-11. 9 502 5976 001 01 Widex Moment Sheer Sound Assist M65

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Android, Google Play and Google Play logo are trademarks of Google Inc.

Subject to change with no further notice.







### WIDEX SOUND ASSIST



SOUND LIKE NO OTHER



## MAKES IT EVEN EASIER TO HEAR THE SOUNDS THAT MATTER TO YOU

Widex Sound Assist is a 'Swiss Army knife' for your hearing aids that fits in the palm of your hand. It can be a table microphone, used as a remote microphone for one-to-one conversations, handle hands-free calls, stream from digital and analogue source, and even be a remote control to adjust your hearing aid settings. Widex Sound Assist also offers Telecoil mode, giving the user the benefit of streaming telecoil audio.



#### Help at your fingertips via the MOMENT app

Via the app it is easy to adjust, for example, volume and sound mixer for the different modes. There are also guides on 'How to pair', 'How to turn it on', and an FAQ section.



#### **TECHNICAL DATA**

- Battery runtime: 10 hrs (mixed use)
- Full charge: 3 hrs
- Li-ion battery
- Micro USB power supply
- Small size: 43x43x16 mm
- Lightweight: 32 grams

