

High risk prevention & intensive treatment of wounds up to grade 4.

When bedsores require intensive treatment, only the best will do. The BariCell system offers effective solutions for a wide range of needs in resident/patient management. The mattress offers multiple configurations to offer the best possible pressure redistribution. The system includes an automatic pump, with alternating, continious low pressure, autofirm and seating function. For bariatric use, we advise the BariCell range for heavy weights from 200 to 450kg.

Activating the Autofirm Mode provides a stable support surface for nursing procedures or to facilitate transfers in or out of bed. The Lock function improves safety and audio-visual Alarm guards pressure levels of the system.



Safety Standards

Complies with: Regulation (EU) 2017 / 745 for medical devices

Protection against Electric shock: Class II , type BF

Protection of Ingress of Liquids:

Mode of Operation: Continuous

Product Specifications

Mattress: Replacement Pump: 230V / 50Hz Width: 100, 120 & 140cm Noise: 0 - 26dB

200 or 220cm Length: Cycle: 10 / 15 / 20 / 25min.

Heigth: 25cm

Design, Quality & Service powered

Patient weight: 450kg by: Medical Solutions